

194mm

103mm

CE IVD REF

SDK60102

Viral Nucleic Acid Isolation Kit

(Silica-Based Spin Column)

30°C

4°C

Storage Temperature: 4~30°C

-10~40°C

-10~40°C

Transportation Temperature: -10~40°C

LOT

Batch code

Date of manufacture

(YYYY/MM/DD)

Use-by date

(YYYY/MM/DD)

CE IVD REF

SDK60102

EC REP

MedNet EC-REP GmbH
Borkstrasse 10 · 48163 Muenster · Germany

bioPerfectus technologies

Jiangsu Bioperfectus Technologies Co., Ltd.

Manufacturer Add: 3rd and 4th floors of Building A(G19), 4th floor of Building F(G14), Ground floor of Building G20, Shuaiyu Village, Fuye village, Sixiang town, Taizhou National Medical, Hi-tech Development Zone, 225300 Taizhou, Jiangsu, PEOPLE'S REPUBLIC OF CHINA.

Tel: 86-21-34637616

Postal Code: 225300

Web: www.bioperfectus.com

Viral Nucleic Acid Isolation Kit

(Silica-Based Spin Column)

Jiangsu Bioperfectus Technologies Co.,Ltd.

50T

Only use for IVD, other details refer to instructions

bioPerfectus technologies

Jiangsu Bioperfectus Technologies Co., Ltd.

Manufacturer Add: 3rd and 4th floors of Building A(G19), 4th floor of Building F(G14), Ground floor of Building G20, Shuaiyu Village, Fuye village, Sixiang town, Taizhou National Medical, Hi-tech Development Zone, 225300 Taizhou, Jiangsu, PEOPLE'S REPUBLIC OF CHINA.

Tel: 86-21-34637616

Postal Code: 225300

Web: www.bioperfectus.com

名称:	病毒核酸提取试剂盒（膜吸附法）外包装英文版		
规格:	50人份/盒		
编号:	TG-000192-00		
版本:	4.0		
尺寸信息:	194*130*103mm（长宽高）		
材质要求:	白卡纸		
印刷要求:	四色彩色印刷		
覆膜要求:	哑膜		
色标要求:	C:80% M:10% Y:0% K:0% C:100% M:0% Y:100% K:0%		
	C:55% M:0% Y:100% K:0% C:100% M:0% Y:100% K:20% C:100% M:0% Y:100% K:30%		
纸质g重:	400G		
存放温度:	室温保存		
备注:	/		
拟制人:	石晓亮		
复核人:	冯永胜		
批准人:	王河清		