easyFLOW TS

transfusion set

features

- 2-way sharp closure-piercing device of the drip chamber with anti-slip surface that guarantees tight and reliable connection with containers/bags
 air inlet with anti-bacterial filter and a coloured (red) flap

- pliable drip chamber
 dropper: 20 drops = 1 ml +/- 0,1 ml
 special blood filter with large area, pore size 200 µm

- soft, flexible tubing
 universal Luer-Lock fitting at the end of tubing
 precise safe flow regulator with integrated holder for the closure-piercing device and with hook for securing the tubing
 name of the manufacturer placed on the flow regulatorboth tips of the device additionally secured with protective caps

- single-use
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil with a coloured code (red) identifying the device

indication

• device used for the administration of blood and blood products

REF	version	length
TS*	with phthalates	150 cm
TS-BF**	without phthalates	150 cm













easyFLOW LINE

extension tube for infusion pump

features

- \bullet tubing made of PVC
- device available in a transparent and amber version for photosensitive medicines
- diameter of the tubing outer: 2,4 mm; inner: 1,24 mm
- New version of the tubing outer : 2,8 mm; inner: 1,8 mm Luer-Lock fitting male and female
- single-use
- non-pyrogenic, non-toxic
- packaging: 1 pc./paper-foil

indication

· device used to extend infusion line

	REF	colour of tube	length	intermediate packaging		
	P450-BF	transparent	45 cm	50 pcs		
	P1500-BF	transparent	150 cm	50 pcs		
	1 1300-DI	transparent	150 CIII	30 pc3		
	PB1500-BF	amber	150 cm	FO		
	PB 1500-BF	amber	150 CIII	50 pcs		
	P2000-BF		200 cm	F.O		
	P2000-BF	transparent	200 Cm	50 pcs		
	PB2000-BF	amber	200 cm	50 pcs		
	PB2000-BF	arriber	200 CIII	50 pcs		



bulk packaging: 24 x 50 pcs

New version of	New version of the tubing outer					
REF	colour of tube	length	intermediate packaging			
P450-18-BF	transparent	45 cm	50 pcs			
D1000 10 DE		100	FO			



