UNIVERSAL SUTURES "HEALING BEYOND COMFORT"





AN ISO & CE CERTIFIED SURGICAL SUTURES AND MESH MANUFACTURING COMPANY We Manufacture complete range of Surgical Sutures, Hernia Mesh & Bonewax

Needle Shape



Needle Type:



International Coding Pattern										
Code 1	Product	Code 2	USP Size	Code 3	Code 4	Needle Curvature	Code 5	Needle Type	Code 6	
1	Braided & Coated Polyglycolic Acid Undyed	A	3		С	3/8 Circle	Z	Round Body / Taper Point		
2	Braided & Coated Polyglycolic Acid	В	2	Ε	D	1/2 Circle	Y	Cutting	F	
3	Monofilament Polyamide (Nylon)	С	1	0	A	Straight	X	Reverse Cutting	E E	
4	Catgut Chromic	D	0	i i	E	5/8 Circle	к	Taper Cutting	⊆	
5	Black Braided Silk	E	2-0	NGTH			т	Trocar Point	LENGTH	
6	Braided Polyester	F	3-0	Ž			R	Blunt Point	U Z	
7	Catgut Plain	G	4-0	<u>Ц</u>			S	Spatula	<u> </u>	
8	Monofilament Polypropylene	н	5-0	SUTURE				•	Ľ ا	
9	Monofilament Polydioxanone	I	6-0	Ē					NEEDLE	
10	Braided & Coated Polyglactin 910	J	7-0	ר					Ž	
11	Monofilament Polyglecaprone 25	к	8-0	1						
12	Braided & Coated Polyglactin 910 Undyed			1						

Example:

 Braided & Coated Polyglycolic Acid
 Size 1
 75cm
 1/2 Circle
 Round Body
 40mm

2 Č 75 Ď Z 40

UNISIL Black Braided Silk Non Absorbable Surgical Suture U.S.P.



UNISIL – Black Braided Silk Suture is a non absorbable, sterile, surgical suture composed of an organic protein called fibroin. This protein is derived from the domesticated species Bombyx mori (B. Mori) of the family Bombycidae.

UNISIL sutures are processed to remove the natural waxes and gums. UNISIL suture is dyed black and coated with a special wax mixture. UNISIL Black Braided Silk is also available in its natural color.

UNISIL Virgin silk is available in which the sericin gum is not removed and serves to hold the filaments together. In-house needle manufacturing provides you the better needles with optimum sharpness which finally delivers smooth riding on tissues during suturing.

Indications:

Black Braided Silk Suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

- Highest purified protein from Species B. Mori results no tissue reaction
- Beeswax coated ultra smooth surface decrease suturing timing followed by no tissue trauma
- Unmatched Knot strength provides satisfaction to user
- USP approved dye increases visibility followed by no reaction


Best Used By:

- Orthopaedic Surgeons,
- General Surgeons,
- Gynaecologists,
- Cardiac Surgeons,
- Endoscopic Surgeons
- **Ultimate Suture For:**
- Skin Closure
- Hard Tissues approximation
- Tendon Repair
- Muscles approximation

Additional Information:

Suture Characteristics	:	Non Absorbable Surgical Suture
Туре	:	Braided
Material	:	Braided Silk
Color	:	Black
Absorption	:	Non Absorbable
USP Range	:	5/0-2
Sterilization	:	EO (Ethylene Oxide)
Shelf Life	:	5 Years
Product Characteristics	:	Excellent tensile strength, Excellent knot tying, Minimum tissue reaction.
Packaging	:	In food grade folder packed in a medical grade pouch, with or without needle
Needle Type	:	Round bodied, cutting edge, taper point, straight, blunt point
MFG. LIC No	:	MFG/MD/2020/000109
OEM. LIC No	:	N. Code: KA/DEVICE/MFG/MD/2020/000109

UNISIL Black Braided Silk Non Absorbable Surgical Suture U.S.P.





U



In Bratislava, on Novem

TC-118 SEC Inte

ber 7th, 2016

Dr. Katarina Srdošová







REGISTERED OFFICE : UNISUR LIFECARE PVT. LTD. No. 15/1,2,3 Andhrahalli Main Road, Acharya Indl. Complex, Vishwaneedam Post, Near Anupama School, BANGALORE - 560091 Karnataka, INDIA Tel.: +91 9108 990 400 | E-Mail: info@universalsutures.com | Web.: www.universalsutures.com

CORPORATE OFFICE :

UNISUR LIFECARE PVT. LTD. Unit No. 303, 3rd Floor, Brigade Rubix, Plot No.: MYS357 in Peenya Plantation, HMT Factory Main Road, BANGALORE - 560013 Karnataka, INDIA Tel.: +91 9108 990 400 | E-Mail: care@universalsutures.com | Web.: www.universalsutures.com

