




BESCO



Hospital Beds



HB036

Manual Homecare Bed

Function:

Back-rest tilting is 0~70°. Knee-rest tilting is 0~45°. Height adjustment operated by crank is 450~760mm.

Standard accessory:

2 pcs of 3" total lock brake castors, 2 pcs of 3" swivel castors. Rectangular tube bed plate, lifting pole base, self-protected metal crank system.

Optional:

Side rails, lifting pole.

Dimension	L2240*W900*H400~580mm
Material	Steel framework
Backrest adjustment angle	0-70°
Leg rest adjustment angle	0-45°
Height adjustment range	400~580mm



HB038

Manual Bed(3-Crank)

Function:

Back-rest tilting is 0~80°. Knee-rest tilting is 0~35°. Height adjustment is 450~730mm.

Material:

Detachable ABS bed board
cold-rolled steel bed plate
self-protected crank system.

Standard accessory:

PE plastic side rails, 5" castor

Optional: Two sectional IV pole, mattress, overbed table, central lock system

Dimension	L2180*W900/990*H450~730mm
Material	Steel framework
Backrest adjustment angle	0-80°
Leg rest adjustment angle	0-35°
Height adjustment range	450~730mm



HB039

Manual Bed

Function:

Back-rest tilting is 0~75°. Knee-rest tilting is 0~45°. Height adjustment is 450~750mm.

Material:

Detachable ABS bed board, cold-rolled steel bed plate, self-protected crank system.

Standard accessory:

Aluminum alloy side rails, 5" castors

Optional: Two sectional IV pole, mattress, overbed table, central lock system

Dimension	L2130*W960*H450~750mm
Material	Steel framework
Backrest adjustment angle	0-75°
Leg rest adjustment angle	0-45°
Height adjustment range	450~750mm