for pulinonary ventilation's devices, modsky CPAP and BILEVEL assisted and signerised.

- Innovative design that guarantees an optimal air-leak control and a high patient comfort level.
 Lightweight.
 Soft replaceable cushion minimizes skin breakdown.
 The kit has all the accessories for positioning.



9

9 1 ı

0



- kit contents:

 * Transparent PVC mask with inflatable and replaceable cushion with inflation valve.

 DEHP free.

 Thumbscrew adjustable forehead cushion with two grips.

 Prote for 0, and CO₂ pressure monitoring and a port for introducing pastrio tableg and with the possibility of introducing a bronchoscope.

 Two connectors 22M and 22F.

 Pump for inflating and deflating the cushion.

 Headgeer with a rapid removal system.

 The mask's connector can be connected to expiratory system (KM-503 not in the kit) provided with a hole for the CO₂ discharge and with a port for the enrichment of O₂.

 Every kit is delivered in a box with a label indicating the contents.

KM-3101	cushion replacement	extra large	KM-310	50
KM-3111	cushion replacement	large	KM-311	50
IOM-3121	cushion replacement	medium	KM-312	50
KM-3131	cushion replacement	small	KM-313	50
KM-3141	cushion replacement	pediatric	KM-314	50



KOO MEDICAL EQUIPMENT

Well, www. Identification.

Well and Identifica