




Urology **Disposables**

Advin Urology
www.advinurology.com



Urology Disposables

Double J Stent / DJ stent / JJ Stent / Ureteric stent / Ureteral Stent / Double Pigtail Ureteral Stent



Double J Stent

Application :

Double J Stent is placed inside the ureter during surgery to ensure drainage of urine from the kidney into the bladder.

Double J stent is made from tecoflex polyurethane which is highly radio opaque and gives better visualization.

Features :

Pigtail shape prevents stent migration.

Superior bio-compatibility.

Thin wall designed for maximum drainage.

Tapered at both the ends - for easy insertion.

Minimized tissue irritation with easy but rapid Traumatic insertion.

Can be used for long term indwelling.

Size	Length	Use Guide	Wire Type	Color
3.0 Fr	16 cms - 26 cms	0.018"	Open End/ Closed End	White/Blue
3.5 Fr	16 cms - 26 cms	0.018"	Open End/ Closed End	White/Blue
4.0 Fr	16 cms - 26 cms	0.025"	Open End/ Closed End	White/Blue
4.5 Fr	16 cms - 26 cms	0.025"	Open End/ Closed End	White/Blue
4.8 Fr	16 cms - 26 cms	0.032"	Open End/ Closed End	White/Blue
5.0 Fr	16 cms - 26 cms	0.035"	Open End/ Closed End	White/Blue
5.5 Fr	16 cms - 26 cms	0.035"	Open End/ Closed End	White/Blue
6.0 Fr	16 cms - 26 cms	0.038"	Open End/ Closed End	White/Blue
7.0 Fr	16 cms - 26 cms	0.038"	Open End/ Closed End	White/Blue
7.0 Fr	16 cms - 26 cms	0.038"	Open End/ Closed End	White/Blue

Types of Loop :

One End multi Loop

Both End Multi Loop

Packing Standard :

Tyvek Pouch as per EN-868-5 series standard.

Different types of Packings :

Double J Stent


Double J Stent with Pusher

Double J Stent with Pusher and Suture


Double J Stent SET (Double J Stent with Pusher, Suture and Guide Wire)

Special Size and length of Double J stents can be available as per buyers request.


!!!!!!



Double J Stent with Pusher



Double J Stent with Pusher and Suture



Double J Stent SET