

Intra Vitreal Injection Set



VKIVK100

Standard Intravitreal Injection Kit

- 1 Eye Speculum
- 1 40 x 40 Drape with Pouch
- 1 70 x 80 Table Drape
- 5 Cotton Swabs
- 4 Gauze Bandages
- 1 30g Cannula
- 1 Scleral Marker
- 1 20g Cannula
- 1 Pair Gloves
- 1 Insulin Injector

VKIVK100

Intravitreal Injection Kit

VKIVK110

Intravitreal Injection Kit

Contents can be determined based upon any desired requirements. Contact us for more information.



VKIVK110

Standard Intravitreal Injection Kit

- 1 Adjustable Eye Speculum
- 1 40 x 40 Drape with Pouch
- 1 70 x 80 Table Drape
- 5 Cotton Swabs
- 4 Gauze Bandages
- 1 30g Cannula
- 1 Scleral Marker
- 1 20g Cannula
- 1 Pair Gloves
- 1 Insulin Injector

Table Cloth

Size can be determined with any desired requirement.



Dimensions (cm)	SMS	SMS Barried	White Nylon	Blue Nylon	Nylon Barried
10.5 x 17.5	-	-	MT10K	-	-
30 x 30	-	-	MT30K	-	-
60 x 60	-	-	MT60K	-	-
75 x 75	-	-	-	MT75N	MT75NB
100 x 100	MT100SMS	MT100SMSB	-	MT100N	MT100NB
100 x 150	MT150SMS	MT150SMSB	-	MT150N	MT150NB
210 x 70	-	-	-	MT210N	MT210NB

Pouched Drape

Dimensions (cm)	One Pouched	Two Pouched
40 x 40	MT100	-
60 x 60	MT102	-
100 x 100	MT10X10	MT10X10D
100 x 130	MT10X13	MT10X13D
100 x 150	MT10X15	MT10X15D
100 x 160	MT10x16	MT10X16D
120 x 130	MT12X13	MT12X13D
120 x 140	MT12X14	MT12X14D
120 x 150	MT12X15	MT12X15D
120 x 160	MT12X16	MT12X16D
130 x 130	MT13X13	MT13X13D
130 x 140	MT13X14	MT13X14D
130 x 150	MT13X15	MT13X15D
130 x 160	MT13X16	MT13X16D
140 x 130	MT14X13	MT14X13D
140 x 140	MT14X14	MT14X14D
140 x 150	MT14X15	MT14X15D
140 x 160	MT14X16	MT14X16D
150 x 130	MT15X13	MT15X13D
150 x 140	MT15X14	MT15X14D
150 x 150	MT15X15	MT15X15D
150 x 160	MT15X16	MT15X16D
150 x 200	MT15X20	MT15X20D
150 x 240	MT15X24	MT15X24D
160 x 130	MT16X13	MT16X13D
160 x 140	MT16X14	MT16X14D
160 x 150	MT16X15	MT16X15D
160 x 160	MT16X16	MT16X16D
160 x 200	MT16X20	MT16X20D
160 x 240	MT16X24	MT16X24D
160 x 250	MT16X25	MT16X25D



Tissue made from SMS. Size can be determined with any desired requirement.


Microscope Drapes




It can be produced suitable for Leica microscope.

MTL500


Microscope Drape




1. Unfold drape, lens cover facing down and away from you, label up.



2. With hands placed inside center fold form pouch opening.




3. Feed drape over furthest protruding microscope components, continue unfolding drape over balance of microscope equipment.



4. Extend balance of drape over microscope arm, gather up extra width with straps.

Ophthalmic Set



MTO100 Ophthalmic Set

Ophthalmic Set contents

- 1 Pouched Drape (non woven)
- 1 Surgical Table Drape (nylon)
- 1 Plain Drape (non woven)
- 2 Gowns
- 2 Towels

Ophthalmic set can be determined based upon any desired requirements.
Contact us for more information.

Drape for Refractive Surgery



Dimensions (cm)	One Pouched	Two Pouched
60 x 60	MTL100	MTL100D
60 x 120	MTL200	MTL200D

MTLR100 Laser Drape Right Eye

MTLL200 Laser Drape Left Eye

Drape made from special nylon especially designed for laser operations.

Operation Gown



MTB100 Operation Gown

MTB101 Operation Gown and Towels