

NON-INVASIVE VENTILATION: C.P.A.P.

CPAP BOUSSIGNAC AND NEBULISATION

Device for respiratory assistance for patients breathing spontaneously allowing nebulisation of broncho-dilators together with a positive pressure.

New

Bachure


▶ **The use of nebulized sympathomimetics** is recommended during acute asthma exacerbations. Indeed this route is more effective and entails fewer side effects than when given intravenously.

▶ **The additive effect of CPAP** allows:

- More distal delivery of the medication
- Reduction of pulmonary resistance
- Decreased work of breathing
- Reduction or even elimination of the Auto-PEEP
- Increased deposition of the nebulized medication by 30%
- Improved administration of broncho-dilators
- Prompt clinical and gasometric improvement.

▶ **Main indication for CPAP and nebulisation: acute asthma**

Make sure to always combine CPAP and conventional medication (beta 2 mimetics and anticholinergic + systemic therapy with corticoids).



CPAP and nebulisation is recommended in case of acute asthma exacerbations (expiratory peakflow < 50%). This technique must be used immediately when the diagnosis is made and the treatment with nebulisation is decided.

Pressure to be applied in the case of acute asthma is 2 to 4 cm H₂O.

The flowrate to be applied on the nebulizer is 6 liters/min.

Renew nebulisation based on clinical improvement of expiratory peakflow and the SpO₂: nebulisation can be renewed every 20 minutes based on the clinical state of the patient.

5577.013 CPAP kit and nebulizer (without mask)

1/24

Contents:

- 1 CPAP device and manometer connector
- 1 nebulizer.

5571 CPAP kit and nebulizer (with mask)

1/24

Contents:

- 1 CPAP device and manometer connector
 - 1 nebulizer
 - 1 face mask (4 sizes)
- code 5571.303 (adolescent)
code 5571.403 (adult - small)
code 5571.503 (adult - medium)
code 5571.603 (adult - large).

5572 CPAP kit and nebulizer (with mask and fixation harness)

1/24

Contents:

- 1 CPAP device and manometer connector
 - 1 nebulizer
 - 1 silicone fixation harness
 - 1 face mask (4 sizes)
- code 5572.303 (adolescent)
code 5572.403 (adult - small)
code 5572.503 (adult - medium)
code 5572.603 (adult - large).

