

Reliable Brands, Innovative Solutions





The gel cuts the air between the skin and the probe, ensuring the transmission of the ultrasonography waves to the device screen more clearly and uninterruptedly In all kinds of ultrasonography and doppler applications. Thanks to its active substances, its high viscosity, and the vacuum technology used during its production, it does not cause any damage the probe, is not affected by the salt in the body during use, and provides a comfortable ultrasonography. Hypoallergenic. Non-irritating to the skin. Free of oil and oily substances, and free of formaldehyde and salt. Not toxic. Odorless.

**FDA approved.** Tests conducted at the University of Michigan have proven to have a better imaging performance than the world's best-known ultrasound gel.

Application Area: Suitable for use in all sonographic device probes. PARABEN FREE





EUROPE'S LEADING MEDICAL GEL MANUFACTURER, THE CHOICE OF MEDICAL PROFESSIONALS!



info@turkuazsaglik.com.tr / www.turkuazsaglik.com.tr