





Arterial



Venous

### Arterial Tubing(red)

| NO. | Description                       | Qty | Mark |
|-----|-----------------------------------|-----|------|
| 1   | Cap for Dialyzer Connector        | 1   |      |
| 2   | Dialyzer Connector(red)           | 1   |      |
| 3   | Cap for Female Luer Lock          | 2   |      |
|     | Female Luer Lock $\phi 4.1$       |     |      |
|     | Clamp(red small)                  |     |      |
|     | Branch Tubing: 2.2*4.1*100mm      |     | A    |
| 4   | Cover for Drip Chamber            | 1   |      |
|     | Dropper                           |     |      |
|     | Drip Chamber (OD:26-30mm,LD:83mm) |     |      |
| 5   | Pump Segment Connector $\phi 4.1$ | 1   |      |
|     | Pump Segment Connector $\phi 5.5$ |     |      |
|     | Pump Segment: 8*12*350mm          |     | D    |
| 6   | Clamp(red small)                  | 1   |      |
|     | Female Luer Lock $\phi 2.5$       |     |      |
|     | Cap for Female Luer Lock          |     |      |
|     | Heparin Line: 1.0*2.5*600mm       |     | B    |
|     | Branch Tubing: 2.2*4.1*50mm       |     | C    |
| 7   | Clamp(red small)                  | 1   |      |
|     | Cap for Female Luer Lock          |     |      |
|     | Transducer Protector              |     |      |
|     | Branch Tubing: 3.5*5.5*500mm      |     | E    |
| 8   | Access Port(red)                  | 1   |      |
| 9   | Clamp(red big)                    | 1   |      |
| 10  | Patient Connector                 | 1   |      |
| 11  | Suspended Spike                   | 1   |      |
| 12  | Main Tubing: 4.5*6.8*750mm        | 1   |      |
| 13  | Main Tubing: 4.5*6.8*600mm        | 1   |      |
| 14  | Main Tubing: 4.5*6.8*200mm        | 1   |      |
| 15  | Main Tubing: 4.5*6.8*1600mm       | 1   |      |

### Venous Tubing(blue)

| NO. | Description                          | Qty | Mark |
|-----|--------------------------------------|-----|------|
| 1   | Cap for Dialyzer Connector           | 1   |      |
| 2   | Dialyzer Connector(blue)             | 1   |      |
| 3   | Cap for Female Luer Lock             | 1   |      |
|     | Transducer Protector                 |     |      |
|     | Clamp(blue small)(blue)              |     |      |
|     | Branch Tubing: 3.5*5.5*500mm         |     | B    |
| 4   | Cover for Drip Chamber               | 1   |      |
|     | Drip Chamber (OD:21.5-30mm,LD:130mm) |     |      |
|     | Filter                               |     |      |
| 5   | Access Port(blue)                    | 1   |      |
| 6   | Cap for Female Luer Lock             | 2   |      |
|     | Female Luer Lock $\phi 4.1$          |     |      |
|     | Clamp(blue small)(blue)              |     |      |
|     | Branch Tubing: 2.2*4.1*100mm         |     | A    |
| 7   | Clamp(blue big)                      | 1   |      |
| 8   | Patient Connector                    | 1   |      |
| 9   | Drainage bag                         | 1   |      |
| 10  | Main Tubing: 4.5*6.8*400mm           | 1   |      |
| 11  | Main Tubing: 4.5*6.8*500mm           | 1   |      |
| 12  | Main Tubing: 4.5*6.8*1800mm          | 1   |      |
|     | Recirculating Connector              | 1   |      |

|             |  |               |  |                                    |  |               |  |                              |   |      |
|-------------|--|---------------|--|------------------------------------|--|---------------|--|------------------------------|---|------|
|             |  |               |  | Prepared by :                      |  | Reviewed by : |  | Description:                 |   |      |
|             |  |               |  | Name                               |  |               |  | Tubing Sets for Hemodialysis |   |      |
|             |  |               |  | Date                               |  |               |  |                              |   |      |
|             |  |               |  | Approved by QA:                    |  |               |  |                              |   |      |
|             |  |               |  | Name                               |  |               |  | CODE:                        |   |      |
|             |  |               |  | Date                               |  |               |  |                              |   |      |
|             |  |               |  | Compatible: Fresenius<br>2008/4008 |  |               |  |                              |   |      |
|             |  |               |  |                                    |  |               |  | BAIN-BL-017E                 |   |      |
|             |  |               |  |                                    |  |               |  | Blood pathway                |   |      |
| CARTON SIZE |  | 530*381*249mm |  | QTY/CTN                            |  | 24sets        |  | Sterilization                | R | DEHP |