


Cystoscopy

Ureterorenoscopes



*High Definition 30.000 PIXEL
Autoclavable, max.134°C/2,2bar/5min*

Ordering n°		BICHLER Ø 4,5/6,5 Fr
46315JA2	Semirigid Ureterorenoscope, 10.000 PIXEL, with 30° angled ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 4,0 Ch for 1x3ch.	
46430JA2	Semirigid Ureterorenoscope, 10.000 PIXEL, with 30° angled ocular, Ø steps by 4,5/6,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 4,0 Ch for 1x3ch.	
Ordering n°		BICHLER Ø 6,0/7,5 Fr
46315RA2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.	
46430RA2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.	
46315RB2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.	
46430RB2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.	
Ordering n°		BICHLER Ø 8,0/9,8 Fr
46315TB2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.	
46430TB2	Semirigid Ureterorenoscope, 30.000 PIXEL, with 30° angled ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.	
Ordering n°		BICHLER Ø 8,5/11,5 Fr
46315VB2	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch.	
46430VB2	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,5/11,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x6 or 2X4,0ch.	
Ordering n°		STRAIGHT Ø 6,0/7,5 Fr
46315RA3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.	
46430RA3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 5°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.	
46315RB3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 5,1 Ch for 1x4 or 2X2,4ch.	
46430RB3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 6,0/7,5Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 5,2 Ch for 1x4 or 2x2,4ch.	
Ordering n°		STRAIGHT Ø 8,0/9,8 Fr
46315TB3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 315mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.	
46430TB3	Semirigid Ureterorenoscope, 30.000 PIXEL, STRAIGHT ocular, Ø steps by 8,0/9,8Ch, direction of view 12°, field of view 90°, WL 430mm, working channel Ø 6,8 Ch for 1x5 or 2X3,0ch.	

ENDOMED