

HITECARE
Hitec Solutions

Epidural kit

Type: AS-E

- Epidural needle 18G(OD: 1.2mm), effective length:80mm, overall length:108mm
The needle tube is marked with distance bands which are 10mm apart and originate from the hub.
- Epidural catheter, OD: 0.8mm, ≥ 900 mm, RADIO-OPAQUE
The catheter is marked with graduation at 10mm intervals from 10mm to 250mm from the tip of patient end. Two bands indicated 100mm, three bands 150mm, four bands 200mm, five bands 250mm.
- Introducer guide
- Catheter connector
- Liquid filter 0.22 μ m
- Lor syringe 10ml

Caution:

Never pull back the catheter through the epidural needle
Do not use if the package is opened or damaged.

BS 6196:1989

Hitec Medical Co., Ltd.

Add: No. 703, Hengnan RD 1328, Minhang District,
201114 Shanghai, PEOPLE'S REPUBLIC OF CHINA



STERILE EO

STERILE



CE 0123



Zhejiang Runqiang Medical Instruments Co., Ltd
599 Rui Feng Street, Gaozhao District, Xiuzhou District, 314031 Jiaxing City,
Zhejiang Province, PEOPLE'S REPUBLIC OF CHINA

EC REP

Shanghai International Holding Corp. GmbH(Europe)
Eiffelstrasse 80, 20537 Hamburg Germany
Tel: 0049-40-2513175 Fax: 0049-40-255726

REF 40103B

LOT 20221209

20221209

20271208

