

Cooperated with Partners from 101 Countries



BSS MEDICAL SUPPLY CO., LIMITED

🌐 No 18, Shabian Road, Torch Hi-tech Industry Zone, Zhongshan, Guangdong, China 528437

☎ +86-760-85280100

✉ info@bssmedical.com

💻 www.bssmedical.com



Wechat



Web



BSS MEDICAL SUPPLY CO., LIMITED

Better Signal Solution Partner

Why BSS?


- **Qifeng Zhou**, invented the conductive gel technology for BSS;
- FDA, CE, ISO13485, Japanese Registration, RoHS, FSC in EU, SFDA, and dozens of patents;
- The daily capacity is 1.5 million pcs by the advanced intelligent automatic production lines;
- Partners from 101 countries, and products quality approved by thousands of international hospitals;
- Full products combination meets variety testing purpose: monitoring, stress test, treadmill, holter, MRI, ICG, ABR;

Professor Qifeng Zhou

Academician of Chinese Academy of Sciences (1999)

The chairman of IUPAC (2018-2020)

President of Peking University (2008-2013)



ECG Electrode

CE ISO FDA SFDA RoHS Pmda JAPAN



Liquid Gel



Hotter



Resting Test



Stress Test



**Pre-wired
Neonatal**



MRI/ X-ray

Better Signal Solution Partner



BSS MEDICAL

2021



Foam Backing

Medical grade material is impervious to fluids, and aggressive high track medical-grade adhesive is ideal for long-term monitoring. Convenient tab allows for easy removal.

Recommended for: ICU/CCU · Surgery · Emergency · Long-Term Monitoring

Shape: Round

Item#	Size(mm)	Description	Package
B5055	50*55	Adult	50/pk, 2000/ cs
B2	45*50	Adult	50/pk, 2000/ cs



B5055

Lot
20,
22

Shape: Rectangle

Item#	Size(mm)	Description	Package
B4	42*46	Adult	50/pk, 2000/ cs
B3244	32*44	Adult	50/pk, 4000/ cs
B3240	32*40	Adult	50/pk, 4000/ cs



B4

Shape: Oval

Item#	Size(mm)	Description	Package
B4346	43*46	Adult	50/pk, 2000/ cs
B6	42*50	Adult	50/pk, 2000/ cs
B3644	36*44	Adult	50/pk, 2000/ cs



B4346

Shape: Round

Item#	Size(mm)	Description	Package
B15	43*43	Adult	50/pk, 2000/ cs
B4040	40*40	Adult	50/pk, 4000/ cs
B8	30*30	Child	50/pk, 2000/ cs
B2525	25.4*25.4	Neonatal	50/pk, 4000/ cs



B15

Soft Material: Non-woven, Micropore

Manufactured with a non-woven, Micropore breathable material that stretches and conforms to the body to promote patient comfort. Porous material minimizes moisture buildup and provides long-term adhesion.

Recommended for: ICU/CCU · Emergency · Long-Term Monitoring

Shape: Round

Item#	Size(mm)	Description	Package
B10	50*50	Adult/ Cloth	50/pk, 2000/ cs
B5050M	50*50	Adult/ Micropore	50/pk, 2000/ cs



B10

Stress Test/ Holter

Medical grade white polyethylene materials is impervious to fluids. It allows sweat to travel around the electrode, making it ideal for stress testing and diaphoretic patients.

Recommended for: Stress Test/ Holter · Treadmill Test · Emergency · ICU/CCU

Shape: Round

Item#	Size(mm)	Description	Package
B7	50*50	Adult/Enhanced Foam	50/pk, 2000/ cs
B6060M	60*60	Adult/ Micropore	50/pk, 2000/ cs



B7

Radiotranslucent

Carbon-based silver chloride sensor allows for an unimpeded view on X-rays.

Recommended for: Radiology · Emergency · Cath Lab · Surgery

Item#	Size(mm)	Description	Package
B4346R	43*46	Adult/Oval	50/pk, 2000/ cs
B5050R	50*50	Adult/ Round	50/pk, 2000/ cs
B3030R	30*30	Child/ Round	50/pk, 2000/ cs



B4346R



Wet / Liquid Gel

Wet Gel electrodes could fast catch the ECG signal. They give dependable quick tracing with easy application and removal.

Recommended for: ·EMS ·Emergency



B7W

Item#	Size(mm)	Description	Package
B5W	43*46	Adult/Oval	50/pk, 1000/ cs
B7W	50*50	Adult / Round	50/pk, 1000/ cs

ABR Test / EEG

To perform an ABR test, electrodes are placed on the surface of the skin at the patient's head. Good signal quality is achieved by preparation of the skin, improving the conductivity. An electrical potential with very low amplitudes is amplified and recorded, while the ABR system delivers acoustical stimuli to the ears.



B19

Recommended for: ·ABR Test ·EEG Test ·Non-invasive EEG

Item#	Size(mm)	Description	Package
B19	32*40	ABR/ EEG	30/ pk, 1200/cs

Pre-wired Neonatal Electrode

The pre-wired neonatal electrode are consist of a conductive adhesive hydro-gel, a Ag/AgCl plated sensing element, a soft cloth or PE Foam substrate with the size of 22mm and a 60 cm flexible tinned copper or carbon lead-wire terminating in a standard DIN connector.



BF22

Item#	Size(mm)	Description	Package
BF22	22*23	Cloth/ Square	1/ pk, 200/cs
BY22	22*22	Cloth/ Round	1/ pk, 200/cs
FF22	22*23	Foam/ Squire	1/ pk, 200/cs
FY22	22*22	Foam/ Round	1/ pk, 200/cs



Indocyanine green(IG)Test

Noninvasive ICG clearance test for estimating hepatic blood flow During halothane and isoflurane anaesthesia.

TEB (ICG) electrodes are designed to work with Thoracic Electrical Bioimpedance (TEB) devices..

Recommended for:

ICU/CCU ·Emergency ·Cardiology ·Anesthesiology



B4090

Item#	Size(mm)	Description	Package
B4045	40*45	Foam	4/pk,100/cs
B4090	40*90	Foam	4/pk,100/cs



B2334

Resting Tab

Tab is used for Resting test to be attached to alligator clip connectors. It has the advantages of fast application and removal from patient's skin.

Item#	Size(mm)	Description	Package
B2334	23*34mm	Ag/AgCl coated PET	100/ pk 10000/cs
B2034	20*34mm	Ag/AgCl coated PET	100/ pk 10000/cs
B2634	26*34mm	Ag/AgCl coated PET	100/ pk 10000/cs

OEM/ODM/ Private label solution

BSS also manufacture wide variety of OEM accounts on a private label basis. We offer the options of Gel, Snap, Backing, Private Label, Product shape, Pouch and Carton Design Solutions.

