


DAHLHAUSEN®**Nebulizer Kits, funnel/funnel**

Nebulizer Kit for inhalation of medicines and liquids by mask and safety tubing with funnel. The set is for single patient use only and should be disposed of after 24 hours.

Aerosol capacity of nebulizer according to EN 13544-1:2007.

Length (oxygen tube): 210 cm

Set content: Mask with Lips style, nebulizer chamber 6 ml with lateral oxygen connection, safety connecting tube with funnel/funnel connector

Sterile: No

Art.no.	Application area
0100008852	Child
0100008853	Adult

P. J. Dahlhausen & Co. GmbH
Emil-Hoffmann-Str. 53
50996 Köln

Tel.: (02236) 39 13 - 0
Fax: (02236) 39 13 - 109