

# Medispo®

### **BEST IN MEDICAL DISPOSABLE**



### WE MAKE THE BEST WE CARE THE MOST

- www.hbmchina.com  $\bigcirc$
- gaochen@hbmchina.com
- 0 1508-09,No.211,GoldenPortRoad,Shanghai,China



## WHAT IS HBM

### HONEST н BETTER B **MUTUAL** Μ

In HBM-We respect the HONESTY, do BETTER, and share MUTUAL benefits.








### BFE **≥ 95%** Type I: Type II: BFE ≥ 98% Type IIR: BFE ≥ 98%

<Bacterial Filtration Efficiency>





**Three Layers Design** 

- Disposable
- Odorless
- 3 layers design
- Comfort and fits perfectly
- Free of fiber glass
- Hypoallergic
- With nose wire support
- Low resistance for better breathing

### **PACKING INFORMATION**

Packing:	50pc/box, 40box/carton
Net Weight:	7.3kgs/ctn
Gross Weight:	7.9kgs/ctn
Carton Size:	52x40x30cm

- and left ear.
- good seal.





Size: Length 175mm Width folded 95mm Width unfolded 170mm Length of Earloop 170mm

### **INSTRUCTIONS**

1. Put on the mask, hang on the loops to right

2. Using hand to mould the nose piece to the shape of your nose, ensure a close fit and

3. Fit check before using, mask sure no air leaks at the mask edge.

### **DISPOSABLE FACE MASK 50PCS/BOX**

Product Code: **TYPE I: 966-EL** TYPE II: 966-EL II TYPE IIR: 966-EL IIR

## **Medispo**<sup>®</sup>

## Anti-Fog **Disposable Medical** Face Mask Type II

### with Ties



## Type II : BFE ≥ 98%

<Bacterial Filtration Efficiency>







### Fluid Resistance

Tests penetration resistance to synthetic blood at three protection levels-120 mm Hg.



### **Bacterial Filtration Efficiency (BFE)**

Tests filtration ability against an aerosol containing bacteria 3.0 microns in size.



### $\approx$ Differential Pressure ( $\Delta$ P)

Gauges air flow from one side of a mask to the other, assessing breathability.

- · Made of non-paper, non-woven material
- Odorless
- Hypoallergic
- · Head straps, of two layers with the same size for perfect fitting on the face
- Ultrasonic Bonding of the mask with the head straps
- Tear resistant.
- Plastic/Wire Nose clip 12cm, for perfect and stable application of the mask on the nose
- Respiratory Resistance Differential Pressure (Delta P ΔP): <60 (Pa / cm<sup>2</sup>).

### **Dimensions:**

Length 175mm Width folded 95mm Width unfolded 180mm controlled-partially lint free

Layer 1 the colored layer is your first line of defense

Laver 2 protects against fluid while maintaining breathability

### Layer 3 soft.non-abrasive inner layer makes the mask comfortable on your skin




### Anti-Fog **DISPOSABLE MEDICAL FACE MASK** TYPE II WITH TIES

### 50CS/BOX

### **PACKING INFORMATION**

Packing: Net Weigl Gross Wei Carton Siz

### Three Layers Design



### Product Code: 966-TO-AF

	50pc/box, 40box/carton
ht:	8.0kgs/ctn
eight:	9.0kgs/ctn
ze:	52x39x52cm