



Adhesive fabric tape SENSIplast

TD-06-I.2.c-5.5

Revision 3

date: 2021-05-26

amendment of: 2024-03-21

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.**ADDRESS:** ul. Pod Borem 18, 41-808 Zabrze, Poland; T: +48 32 376 07 00-02,06 F: +48 32 370 38 08; www.zarys.com**CLASSIFICATION**

wyrób medyczny klasy I zgodny z Rozporządzeniem PE i Rady (UE) 2017/745 (MDR)

INTENDED PURPOSE

- device intended for stabilisation, protection and fixation of drains, catheters, medical tubes not requiring specialised fixation at their place of application as well as primary dressings (contact), non-adhesive, applied directly to wounds and skin damages

FEATURES

- adhesive on a cotton fabric backing
- white colour
- wound on a spool
- covered with an acrylic adhesive applied by the continuant method
- can be divided lengthwise and crosswise without the use of scissors
- serrated edge finishing
- good air and water vapour permeability
- good adhesion
- thin, flexible, conforms to the body contours, guarantees freedom of movement
- hypoallergenic
- single-use
- non-sterile



CATALOGUE NUMBER	DIMENSIONS	INTERMEDIATE PACKAGING (SALES) BOX	BULK PACKAGING CARTON
STA12505W	1,25 cm x 5 m	24 x 1 pc.	12 x 24 pcs
STA12509W	1,25 cm x 9,14 m	24 x 1 pc.	12 x 24 pcs
STA25005W	2,5 cm x 5 m	12 x 1 pc.	12 x 12 pcs
STA25009W	2,5 cm x 9,14 m	12 x 1 pc.	12 x 12 pcs
STA50005W	5 cm x 5 m	6 x 1 pc.	12 x 6 pcs
STA50009W	5 cm x 9,14 m	6 x 1 pc.	12 x 6 pcs
STA25005W-L*	2,5 cm x 5 m	12 x 1 pc.	12 x 12 pcs

*version with stronger acrylic glue