

Item Code 910004

INTERNATIONAL
PATENT

Volume Accuracy
 $\pm 2,5\%$ or 50 ml



FlowMIR™

Disposable Turbine Flowmeter
MIR exclusive product
International patent

Each turbine, which includes a cardboard mouthpiece, has been individually factory tested with a computerized system

For each patient, after the Spirometry test, both turbine and mouthpiece are thrown away



FlowMIR Dispensers (10 and 60 pcs)



Comfortable
Packaging



Singularly Tested
and Packed



Eliminates staff
clean-up time



No Sterilisation



No Calibration



No Cross
Contamination



No Antibacterial
Filter



Always 100%
accurate and hygienic



ATS/ERS Compliant
for accuracy



Not affected by
Vapour Condensation



Not affected by
Ambient conditions



The Best Sensor
for Spirometry

Consumables / Accessories



Paper roll for printer

Code 910350

MIR Model Spirolab

Standard thermal paper roll for printer.



Plastic noseclip

Code 910320

MIR Model All Spiro

Plastic noseclip.



Calibration syringe 3L

Code 919000

MIR Model All Spiro

3L calibration syringe.



Oximetry finger sensor

Code 919024

MIR Model All devices with oximetry option

Reusable, ADULT.



Oximetry finger sensor

Code 939005

MIR Model All devices with oximetry option

Reusable, PAEDIATRIC 5-45 kg.