## Runibeda

Gamyba: Vilniaus g. 28a, Prienai LT-59115, Lietuva; Biuras: Europos pr. 96, Kaunas LT-46351, Lietuva tel. +370 673 54050, el. p. <u>info@runibeda.com</u>, kodas 305573341, PVM mokėtojo kodas LT100013424612 AB,,SEB", B.k.: 70440, A.s.: LT777044090100399773

> Transportation stretcher "LECTICA" ZG2-030 Basic UDI-ID: (477905279LECTICX4)

- Description Stretcher intended for the transport of patients
- Hydraulic height adjustment, height adjustment on both sides, and both corners with hydraulic pump with a foot pedals.
- Hydraulic stretcher type, 4- Sections.
- Electrostatically painted stainless steel material. Stretcher have smooth surface, can be easy cleaned and disinfected.
- The reclining positions (Trendelenburg, anti-Trendelenburg) are hydraulically adjustable with the foot pedals.
- The positions (height, trendelenburg, antitrendelenburg) are adjusted by foot pedals on each side of the longer side of the device.
- Backrest adjustment assisted by gas spring
- Foam mattress with Removable, anti-static and fireproof, waterproof, and compatible with the size of the stretcher, 90mm mattress thickness. Mattress antistatic cover can be disinfected and washable. Mattress fixed to the frame with Velcro fastenings.
- Transparency on x-rays (x-ray, arm C). The lying surface X-ray permeable and allow an X-ray examination of the patient lying through the lying surface from the lungs, pelvis, and lower extremities. Niches for x-ray cassettes 4pcs..X-raythere with base on the bed for the cassette.
- Equipped with vertical foldable, lifting and lowering side rails. Length of the side supports 2/3 of the total length of the frame. The side supports positioned so that the patient's head remains open (closer to the foot) from the side supports. Side supports lowered below the mattress to ensure unrestricted access to the patient.
- Durable handles at both ends that allows to push the device comfortably.
- The head end of the patient has one height-adjustable drip leg of strong construction.
- Wheel corners are equipped with circular bumpers. Protective corner bumpers/wheels -Shockproof protection on all four corners made of soft plastic material.

## Runibeda

Gamyba: Vilniaus g. 28a, Prienai LT-59115, Lietuva; Biuras: Europos pr. 96, Kaunas LT-46351, Lietuva tel. +370 673 54050, el. p. <u>info@runibeda.com</u>, kodas 305573341, PVM mokėtojo kodas LT100013424612 AB,,SEB", B.k.: 70440, A.s.: LT777044090100399773

- Anti-static and anti-tracing, 360° Swivel, 200mm Castors with central braking mechanism and steering control. System is activated with the foot pedals. Central brake pedals for four wheels. Central locking system has three functions: all wheels locked, steering wheel for direction drive, wheels unlocked.
- The device has a fifth middle wheel. The total of 4 + 1 wheels.
- Each corner has bumpers or guards to protect the frame. There are no parts projecting laterally beyond the side supports (e.g. knobs, handle), except for the bumpers.
- Oxygen cylinder attachment, in the middle of the frame at the bottom. Cylinders suitable for the oxygen cylinder holder with dimensions ranging from: d = 100 to 140 mm, length 420 to 900 mm.
- EC declaration.
- Warranty- 24 months.

## Technical Details

Length : 2220mm Width : 830mm Backrest Angle : 82° adjustable with gas spring TR/RTR Angle : +17°/-17° (with a shock absorber) Knee Angle : 40° Foot End Angle : 30° Height Range(without mattress) 550 – 800mm Castors Diameter : 200 mm Mattress Surface : 75 x 195 cm 5th Castor Diameter : 150mm





Gamyba: Vilniaus g. 28a, Prienai LT-59115, Lietuva; Biuras: Europos pr. 96, Kaunas LT-46351, Lietuva tel. +370 673 54050, el. p. <u>info@runibeda.com</u>, kodas 305573341, PVM mokėtojo kodas LT100013424612 AB"SEB", B.k.: 70440, A.s.: LT777044090100399773

Mattress Thickness : 90mm

Safe Working Load Capacity : 310kg.



