



## Order Information

ORDER CODE	DESCRIPTION	CARTRIDGE COLOR	STAPLE SIZE (INNER TO OUTER ROW)	QUANTITY
EGIA30AVM	Endo GIA™ 30 mm Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2 mm, 2.5 mm, 3 mm	Pack 6
EGIA30AMT	Endo GIA™ 30 mm Articulating Medium/Thick Reload with Tri-Staple™ Technology	Purple	3 mm, 3.5 mm, 4 mm	Pack 6
EGIA30CTAVM	Endo GIA™ 30 mm Curved Tip Articulating Vascular/Medium Reload with Tri-Staple™ Technology	Tan	2 mm, 2.5 mm, 3 mm	Pack 6
EGIA30AV*	Endo GIA™ 30 mm - Gray Reload	Gray	2 mm, 2 mm, 2 mm	Pack 6

## OTHER CHARACTERISTICS & REQUIREMENTS

\* EGIA30AV, while offering a slimmer fixed anvil, is a vascular reload that does not come with Tri-Staple™ technology

## OVERVIEW

### Designed to Navigate and Staple Tissue in Tight Spaces

The 30 mm reload with Tri-Staple™ technology is designed for use in situations where space is limited. Its extra-short profile provides access and maneuverability in the tightest of cavities while still delivering the performance of Tri-Staple™ technology. The 30 mm reload is also available in a curved tip option for access to vessels and anatomical structures in confined spaces.

[Order Information >](#)

#### About Medtronic

[Medtronic Corporate ANZ](#)

[Privacy Statement](#)

[Terms of Use](#)

[Terms & Conditions of Supply](#)

[External Research Program](#)

#### Minimally Invasive Therapies

[Products](#)

[Terms & Conditions of Supply](#)

[Sitemap](#)

[Australia/New Zealand Customer Service](#)



[Contact Us](#)

#### Country Office

Level 3 - Building 5  
Central Corporate Centre  
666 Great South Road  
Penrose, Auckland  
New Zealand