

HANGZHOU SHANYOU MEDICAL EQUIPMENT CO., LTD
No. 138, Louta Development Zone, Guancun Village, Louta
Town,
Xiaoshan District, Hangzhou, Zhejiang 311266, China
Tel: 86-571-82304593 Fax: 86-571-82305195

Standard Endotracheal Tube

The endotracheal tube is inserted into the trachea through glottis to establish and maintain an artificial airway when take the anesthesia and emergency treatment.



Features:

1. Made of clear, soft, medical-grade thermosensitive PVC, high quality
2. Marking line in cm
3. With Murphy eye to avoid complete respiratory obstruction
4. Radiopaque line throughout the tube for X-ray Visualization, Transparent
5. Sterile, EO
6. Individually packaged
7. CE, ISO Certified
8. Well placed connector

Standard Endotracheal Tube	
Uncuffed	
Size	Product No.
2	ET1P20
2.5	ET1P25
3	ET1P30
3.5	ET1P35
4	ET1P40
4.5	ET1P45
5	ET1P50
5.5	ET1P55
6	ET1P60
6.5	ET1P65
7	ET1P70
7.5	ET1P75
8	ET1P80
8.5	ET1P85
9	ET1P90
9.5	ET1P95
10	ET1P10



Hangzhou Shanyou Medical Equipment Co., Ltd.

Editor: Amira Ma (马园晓) 马园晓