



631-1163

SILANE COATED, COLOUR-FROSTED SLIDES

Silane coated (aminoalkylsilane) glass slides producing a sticky surface, which can bind tissue sections to glass.

- Smooth edges for operational safety
- Marking with pencils, marker pens and conventional labelling systems
- Shelf life of 18 months

Recommended for H&E stain, immunohistochemistry (IHC) and hybridisation on frozen sections in histopathology examination.

50 pieces per box. 50 boxes per case.

Fragile items should be stored in dry and proper laboratory temperature conditions. Not tested for use in microarray assays.

Description	WxD (mm)	Thickness (mm)	Colour	Material	Pk	Cat. No.
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Blue	Glass	2.500	631-1163
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Green	Glass	2.500	631-1164
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Orange	Glass	2.500	631-1165
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Pink	Glass	2.500	631-1166
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Yellow	Glass	2.500	631-1167
Silane coated colour-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Lilac	Glass	2.500	631-1168



COLOUR-FROSTED MICROSCOPE SLIDES

Microscope slides recommended for haematoxylin and eosin (H&E) stain, immunohistochemistry (IHC) and hybridisation on frozen sections in histopathology examination.

- Smooth edges for operational safety
- Sticky surface for stronger adhesion than with poly-L-lysine
- Marking with conventional labelling system, pencils, marker pens and thermal transfer printers

50 pieces per box. 50 boxes per case.

To be stored in dry and proper laboratory temperature conditions – shelf life of 60 months.

Description	WxD (mm)	Thickness (mm)	Colour	Material	Pk	Cat. No.
Colour-frosted microscope slides	76x25	1,0 (-0,05/+0,1)	Blue	Glass	2.500	631-1169
Colour-frosted microscope slides	76x25	1,0 (-0,05/+0,1)	Green	Glass	2.500	631-1172
Colour-frosted microscope slides	76x25	1,0 (-0,05/+0,1)	Orange	Glass	2.500	631-1173
Colour-frosted microscope slides	76x25	1,0 (-0,05/+0,1)	Pink	Glass	2.500	631-1175
Colour-frosted microscope slides	76x25	1,0 (-0,05/+0,1)	Yellow	Glass	2.500	631-1176

INKJET PRINTER SLIDES

Microscopy slides are designed with colour-frosted ends for H&E stain on frozen sections in histopathology examination.

- Shelf life of 60 months
- Pre-cleaned, without mark, print or contamination risks
- Marking with conventional labelling system, pencil, marker pens and inkjet printers

Specially produced for all inkjet printers.

50 pieces per box. 50 boxes per case.

Fragile items should be stored in dry and proper laboratory temperature conditions.



Description	LxW (mm)	Thickness (mm)	Colour	Material	Pk	Cat. No.
Inkjet printer coloured-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	White	Glass	2.500	631-1177
Inkjet printer coloured-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Blue	Glass	2.500	631-1178
Inkjet printer coloured-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Green	Glass	2.500	631-1179
Inkjet printer coloured-frosted slides on 20 mm	76x25	1,0 (-0,05/+0,1)	Yellow	Glass	2.500	631-1182