

Accessories

HUMIDIFIERS AND NEBULISERS

APPLICATIONS

For use in insufflation or inhalation.

Humidifier bubble jet for the moistening of e.g. oxygen, compressed air or nebuliser unit for the nebulisation (aerosol treatment) of sterile water or medications.

TECHNICAL DATA

Design:	Chrome-plated all-metal design. Suitable for connection to a flowmeter.
Material:	Head: brass, chrome-plated Bottle: Apec®
Medium:	Sterile water
Inlet:	9/16-18 UNF
Capacity:	250 ml total volume
Outlet:	6 mm flow nozzle (humidifier) or 22 mm ventilation hose connector (nebuliser)
Dimensions (WxHxD) /Weight:	90x190x55 mm / 320 g (humidifier) 90x190x55 mm / 370 g (nebuliser)



HUMIDIFIER / NEBULISER

Humidifier unit, complete	904.836
Nebuliser unit, complete	904.837
Bottle, 250 ml, imprinted	900.922

FLOWMETER SPARE PARTS

Hose connector for flowmeter with 9/16" union nut	900.619
Measuring tube, 0-6 l/min incl. float	326.458
Measuring tube, 0-15 l/min incl. float	800.278
Measuring tube, 0-30 l/min incl. float	800.279
Measuring tube sleeve	100.089
O-ring 18.77x1.78 NBR 70	100.014
O-ring 6 x 2 NBR 70	100.083
Tube Flowmeter repair kit	900.430
Kolibri flowmeter repair kit	900.431
Med. flowmeter servicing fee	902.047