



| Order Information | | | | |
|-------------------|---|-----------------|----------------------------------|----------|
| ORDER CODE | DESCRIPTION | CARTRIDGE COLOR | STAPLE SIZE (INNER TO OUTER ROW) | QUANTITY |
| EGIA45AXT | Endo GIA™ 45 mm Articulating Extra Thick Reload with Tri-Staple™ Technology | Black | 4 mm, 4.5 mm, 5 mm | 1 Each |
| EGIA60AXT | Endo GIA™ 60 mm Articulating Extra Thick Reload with Tri-Staple™ Technology | Black | 4 mm, 4.5 mm, 5 mm | 1 Each |

OVERVIEW

An advanced solution for thick tissue applications

The first reload designed specifically for use in extra-thick tissue applications and the only extra-thick reload available in 45 mm and 60 mm linear lengths, the Endo GIA™ black reload with Tri-Staple™ technology is ideal for thick, dense tissues commonly found in challenging bariatric and thoracic procedures.

[Order Information >](#)

TECHNOLOGY

Advanced technology for optimal interaction between stapler and tissue

No single stapler can address the wide range of tissue management issues surgeons face on a daily basis. That's why we have developed a portfolio of stapling solutions designed to meet specific needs and challenges. Each reload incorporates proprietary features that deliver performance conventional staplers simply can't match.

All of our stapling reloads incorporate our proprietary Tri-Staple™ technology with its design elements including a strong fixed anvil, stepped cartridge face and staggered staple heights.

[LEARN MORE >](#)

Corporate Resources

[Medtronic.co.uk](#)

[Terms Of Use](#)

[Privacy Statement](#)

Medical Surgical Portfolio

[Clinical Solutions](#)

[Clinical Education](#)

[Email Preference](#)

Additional Resources

[Animal Health](#)



[Contact Us](#)

Medtronic Limited

Building 9, Croxley Park
Hatters Lane
Watford, WD18 8WW
United Kingdom