

MONOLAC

MONOLAC is a mid-term absorbable, synthetic, monofilament suture made of glycolide and ϵ -Caprolactone, available in natural colour for the closure of superficial wounds and in violet for all other indications, when a mid-term monofilament is needed.

During absorption the material elicited only a very slight inflammation.

The material is degraded into adipic acid, which is completely metabolised in the body. The process begins with a loss of tensile strength, followed by a loss of mass at a later stage.



UK03
53

Absorbable - Monolac

Indication:	Soft tissue in all specialities, Skin Closure (e.g. intracutaneous, subcutaneous), Plastic, Gynaecology, Urology and Gastrointestinal surgery.						
Contraindication:	Monolac should not be used in tissues subject to expansion, stretching or distension, or which might require long-term mechanical support.						
Sizes:	6/0 to 2 USP						
Structure:	Monofilament						
Absorption method:	Hydrolysis						
Strength Retention Profile:	MID-TERM HEALING (30 days) <table> <tr> <td>1 Day</td><td>100%</td></tr> <tr> <td>7 Days</td><td>60% violet/50% natural</td></tr> <tr> <td>14 Days</td><td>30% violet/20% natural</td></tr> </table>	1 Day	100%	7 Days	60% violet/50% natural	14 Days	30% violet/20% natural
1 Day	100%						
7 Days	60% violet/50% natural						
14 Days	30% violet/20% natural						
Mass Absorption profile:	90 - 120 days						

Lot 14





MONOLAC

Needle	EP	USP	Length	Colour	Cat. No	Other information
--------	----	-----	--------	--------	---------	-------------------

UK03
55

HR: ½ circle, round body needle

2 x HR 30 	3	2/0	1 x 90 cm	Violet	MV 1373	
HR 35 	3 3.5	2/0 0	1 x 75 cm 1 x 75 cm	Violet Violet	MV 0271 MV 0272	
HR 40 	2 4	3/0 1	1 x 90 cm 1 x 90 cm	Violet Violet	MV 0373 MV 0376	
	3.5 4	0 1	1 x 75 cm 1 x 75 cm	Violet Violet	MV 0276 MV 0277	
	2 4 5	3/0 1 2	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	MV 0377 MV 0397 MV 0398	

HR s: ½ circle, round body needle - slim

HR s 18 	1.5	4/0	1 x 75 cm	Violet	MV 1256	Slim needle
HR s 22 	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	MV 1261 MV 1262	Slim needle
HR s 22 Black 	2	3/0	1 x 75 cm	Violet	MV 2262	Slim, black needle
HR s 27 	3	2/0	1 x 75 cm	Violet	MV 1267	Slim needle

Absorbable - Monolac



Label 15

MONOLAC violet ligatures

UK03

54

Pre-cut threads

EP	USP	Length	Cat. No
3	2/0	1 x 150 cm	MV 0085
3.5	0	1 x 150 cm	MV 0086



MONOLAC

Absorbable - Monolac

Needle	EP	USP	Length	Colour	Cat. No	Other information
HR: ½ circle, round body needle						
 HR 13	0.7	6/0	1 x 75 cm	Violet	MV 0251	
	1	5/0	1 x 75 cm	Violet	MV 0252	
	1.5	4/0	1 x 75 cm	Violet	MV 0253	
 HR 17	1.5	4/0	1 x 75 cm	Violet	MV 0254	
 HR 18	1	5/0	1 x 75 cm	Violet	MV 0256	
	1.5	4/0	1 x 75 cm	Violet	MV 0257	
	2	3/0	1 x 75 cm	Violet	MV 0258	
 HR 22	1.5	4/0	1 x 75 cm	Violet	MV 0261	
	2	3/0	1 x 75 cm	Violet	MV 0262	
	3	2/0	1 x 75 cm	Violet	MV 0263	
 2 x HR 22	3	2/0	1 x 45 cm	Violet	MV 0963	
	1.5	4/0	1 x 75 cm	Violet	MV 0965	
	2	3/0	1 x 90 cm	Violet	MV 0964	
 HR 26	2	3/0	1 x 75 cm	Violet	MV 0166	
	3	2/0	1 x 75 cm	Violet	MV 0167	
 HR 27	1	5/0	1 x 75 cm	Violet	MV 0264	
	1.5	4/0	1 x 75 cm	Violet	MV 0265	
	2	3/0	1 x 75 cm	Violet	MV 0266	
	3	2/0	1 x 75 cm	Violet	MV 0267	
 2 x HR 27	2	3/0	1 x 90 cm	Violet	MV 0366	
			1 x 90 cm	Violet	MV 0367	
	2	3/0	1 x 75 cm	Violet	MV 2266	
 HR 30	1.5	4/0	1 x 90 cm	Violet	MV 2275	
	2	3/0	1 x 90 cm	Violet	MV 2276	
 HR 30	2	3/0	1 x 75 cm	Violet	MV 0269	
	3	2/0	1 x 75 cm	Violet	MV 0273	
	3.5	0	1 x 75 cm	Violet	MV 0274	