



PS301-OB-EN
SAM ST PELVIC IMMOBILIZER

81 - 127 cm

Standard 81 to 127 cm

Large 91 to 137 cm

Small 69 to 114 cm

Designed to ensure the correct force on pelvic fractures, the patented "Autostop" system does not allow the compression force to exceed 14.97 kg.

This factor is vital for highly stressful situations in which medical staff could put too much pressure on the patient and be harmful to them.

It provides the ideal force at all times, eliminating any uncertainty about how much to tighten the pelvic belt. The system is programmed so that the buckle stops pressing when the correct compression force is reached.

The two prongs of the buckle come loose and prevent the belt from tightening further.

