
UNIVERSAL SUTURES "HEALING BEYOND COMFORT"



Sutures Mesh Bone Wax

Needle Shape



Needle Type:











| International Coding Pattern | | | | | | | | | | |
|------------------------------|---|--------|----------|--------------|--------|------------------|--------|--------------------------|--------|--|
| Code 1 | Product | Code 2 | USP Size | Code 3 | Code 4 | Needle Curvature | Code 5 | Needle Type | Code 6 | |
| 1 | Braided & Coated Polyglycolic Acid Undyed | Α | 3 | | С | 3/8 Circle | Z | Round Body / Taper Point | | |
| 2 | Braided & Coated Polyglycolic Acid | В | 2 | LENGTH in cm | D | 1/2 Circle | Υ | Cutting | u u | |
| 3 | Monofilament Polyamide (Nylon) | С | 1 | | Α | Straight | Х | Reverse Cutting | | |
| 4 | Catgut Chromic | D | 0 | | E | 5/8 Circle | К | Taper Cutting |] ≗ | |
| 5 | Black Braided Silk | E | 2-0 | | | | T | Trocar Point | | |
| 6 | Braided Polyester | F | 3-0 | | | | R | Blunt Point | LENGTH | |
| 7 | Catgut Plain | G | 4-0 | | | | S | Spatula | | |
| 8 | Monofilament Polypropylene | Н | 5-0 | I RE | | | | |] " | |
| 9 | Monofilament Polydioxanone | I | 6-0 | SUTUR | | | | | NEEDLE | |
| 10 | Braided & Coated Polyglactin 910 | J | 7-0 | S | | | | | 🖁 | |
| 11 | Monofilament Polyglecaprone 25 | К | 8-0 | | | | | | | |
| 12 | Braided & Coated Polyglactin 910 Undyed | | | 1 | | | | | | |

Example:

Braided & Coated Polyglycolic Acid Size 1 75cm 1/2 Circle Round Body 40mm

2 C 75 D Z 40

UNIGLYDE MONO

Monofilament Polyglecaprone 25
Absorbable Surgical Suture (Synthetic)



UNIGLYE MONO - Polyglecaprone 25 suture is a monofilament synthetic absorbable surgical suture prepared from a copolymer of glycolide and epsilon-caprolactone. Monofilament Poliglecaprone 25 suture copolymer has been found to be non antigenic, non pyrogenic and elicits only a slight tissue reaction during absorption.

In-house needle manufacturing provides you the better needles with optimum sharpness which finally delivers smooth riding on tissues during suturing.

Indications:

Poliglecaprone 25 sutures are indicated for use in general soft tissue approximation and/or ligation, but not for use in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

Ultimate Suture For:

- Female's skin approximation
- Subcutaneous approximation

Best Used By:

- · Plastic Surgeons,
- · Gynaecologists,
- · Tricologists,
- · General Surgeon,
- Orthopaedics

Additional Information:

Characteristics : Synthetic Absorbable Surgical Suture

Type : Monofilament

Material : 75%poly (glycolide) – 25% co-caprolactone synthetic absorbable surgical

suture

Coating : -

Color : Undyed / Violet Wound Support : 21 to 28 days

Absorption : Absorption takes place through hydrolysis completely in 90 to 120 days.

Tensile Strength : 50% tensile strength is retained for 14 days.

USP Range : 4/0 - 2

Sterilization : EO (Ethylene Oxide)

Shelf Life : 3 Years

Product Characteristics : Excellent tensile strength, Excellent knot tying, Minimum tissue reaction

Packaging : In aluminium foil packages with or without needle

Needle Type : Round bodied, cutting edge, taper point, straight, blunt point

MFG. LIC No : MFG/MD/2020/000109

OEM. LIC No : N. Code: KA/DEVICE/MFG/MD/2020/000109

UNIGLYDE MONO

Monofilament Polyglecaprone 25
Absorbable Surgical Suture (Synthetic)



| Code No. | SUTURE | NEEDLE DESCRIPTION |
|----------|-----------------|--------------------|
| | USP SIZE LENGTH | |

| Polyglecaprone 25 Undyed Codes | | | | | | | | |
|-------------------------------------|-----|--------|---|--|--|--|--|--|
| S 1614 | 0 | 150 cm | Without Needle | | | | | |
| N 1205 | 4-0 | 70 cm | 3/8 Circle Reverse Cutting 16mm | | | | | |
| N 1326 | 3-0 | 70 cm | 3/8 Circle Cutting 25mm PRIME | | | | | |
| N 1326A | 3-0 | 45 cm | 3/8 Circle Cutting 25mm PRIME | | | | | |
| N 3326 | 3-0 | 70 cm | 3/8 Circle Reverse Cutting 26mm | | | | | |
| N 1641 | 2-0 | 70 cm | 1/2 Circle Round Bodied 30mm | | | | | |
| N 1642 | 1-0 | 70 cm | 1/2 Circle Round Bodied 30mm | | | | | |
| N 1651/35 | 2-0 | 35 cm | Straight Cutting 60mm | | | | | |
| N 3203 | 5-0 | 45 cm | 3/8 Circle Reverse Cutting 13mm | | | | | |
| N 3204 | 5-0 | 45 cm | 3/8 Circle Reverse Cutting 16mm | | | | | |
| Polyglecaprone 25 Dyed Violet Codes | | | | | | | | |
| N 4228M | 2 | 90 cm | 1/2 Circle Round Bodied 45mm | | | | | |
| N 4241M | 2-0 | 76 cm | 1/2 Circle Round Bodied 30mm | | | | | |
| N 4242M | 0 | 76 cm | 1/2 Circle Round Bodied 30mm | | | | | |
| N 4246M | 0 | 76 cm | 1/2 Circle Round Bodied 40mm | | | | | |
| N 4259M | 1 | 76 cm | 1/2 Circle Round Bodied 40mm | | | | | |
| S 2214M | 0 | 152 cm | Without Needle | | | | | |
| S 2215M | 1 | 152 cm | Without Needle | | | | | |
| N 2214M | 0 | 152 cm | 1/2 Circle Round Bodied 30mm | | | | | |
| N 2215M | 1 | 152 cm | 1/2 Circle Round Bodied 40mm | | | | | |
| N 4226M | 1 | 76 cm | 3/8 Circle Round Bodied 45mm | | | | | |
| N 4273M | 1 | 76 cm | 1/2 Circle Reverse Cutting 50mm | | | | | |
| N 4245M | 2-0 | 76 cm | 1/2 Circle Round Bodied 40mm | | | | | |
| N 4248M | 1 | 100 cm | 1/2 Circle Round Bodied 50mm Heavy | | | | | |
| N 4255M | 2-0 | 76 cm | 1/2 Circle Cutting 30mm Heavy | | | | | |
| N 4258M | 0 | 76 cm | 1/2 Circle Round Bodied 40mm Heavy | | | | | |
| N 1326DY | 3-0 | 70 cm | 3/8 Circle Cutting 25mm PRIME | | | | | |
| N 1326A DY | 3-0 | 45 cm | 3/8 Circle Cutting 25mm PRIME | | | | | |
| N 1326A DY | 3-0 | 45 cm | 3/8 Circle Cutting 25mm PRIME | | | | | |
| 11F70EZ26DA Lot 191 | 3-0 | 70 cm | 5/8 Circle Round Body 26cm Double Armed | | | | | |