

Alginate



Alginate Wound Dressing

UMT Alginate Wound Dressing SEQUA™sorb is an advanced wound dressing manufactured from sodium alginate extracted from natural seaweeds. When comes in contact with a wound exudates, calcium in the dressing is exchanged with sodium from wound fluid turning the dressing into a gel which maintain a moist wound healing environment and is ideal for the management of exuding wounds.



LMT

Foshan United Medical Technologies Ltd

Company Address : Southern Medical Devices Industrial Park,
89 Taoyuan East Road, Shishan, Nanhai, Foshan,
Guangdong Province, 528225 China

Tel: +86-757-89959856

Fax: +86-757-89959013

E-mail: info@umtmedical.com **Web:** www.umtmedical.com



Features and Benefits

- High Absorbency
- Strong Wet Strength
- Gel-forming Property
- Pain Reducing
- Cost Reducing
- CE Marked



before

after

Indications:

Alginate can be used in multiple wound types, including diabetic wounds, venous wounds, pressure ulcers, cavity wounds etc.



Product Specification

5X5	10 pcs /carton
10X10	10 pcs /carton
5X15	10 pcs /carton
10X20	10 pcs /carton
20X30	10 pcs /carton
10X12.5	10 pcs /carton
15X15	10 pcs /carton
15X25	10 pcs /carton
2.5X30	10 pcs /carton
2.5X45	10 pcs /carton