

1. PRODUCT

UNIVERSAL DRAPE PACK

2. EMDN CODE

T02010101

Product Code: 200-001-S.1

2. Product conformities:

ISO 9001, ISO 13485, EN 13795, ISO 11135:2014, ISO 14644-3:2019, ISO 14791:2012, ISO 14937: 2009,
ISO 10993-7:2008

3. DESCRIPTION

Mr.Baki Polat Canavar


ZEYNI MEDİKAL .
TEKSTİL İNŞAAT SAN. VE TİC.LTD.ŞTİ.
Yunusemre Mah. 12. Yıldırım Sk. No:17
Yıldırım / BURSA / TÜRKİYE
Tic Sic No: 095814140000001
Yıldırım Mah. 12. Yıldırım Sk. No:17

4. STORAGE

+5°C to + 25°C, in a dry and humid environment, protect from direct sunlight.

5. VALIDITY

3 Years after production

Zeyni Medikal Tekstil Co. Ltd.

Yunusemre Mah. 12. Yıldırım Sk No.17 16270 Yıldırım Bursa Türkiye
Tel :+90224 369 3232 Mob:+90 533 143 5356
E:baki@zeynimedical.com
W:www.zeynimedical.com.tr



Packing Details



Zeyni Medikal Tekstil Co. Ltd.

Yunusemre Mah. 12. Yıldırım Sk No.17 16270 Yıldırım Bursa Türkiye
Tel : +90224 369 3232 Mob: +90 533 143 5356
E: baki@zeynimedical.com
W: www.zeynimedical.com.tr



Product Specifications

Color	Blue
Fabric	Nonwoven
Sterile	Yes
Disposable	Yes
Latex Free	Yes
Skin Friendly	Yes
Recycleable	Yes
Reinforcement Type	Fully Impervious
Alcohol Repellent Material	Level 10
Antiseptic Present	Yes
Without rosin	Yes
Individually packed	Yes



1. Reinforced mayo stand cover 70-90x135-145cm, absorbent zone 60-70x80-90cm with guiding pictograms, waterproof
2. Absorbent adhesive drape 85-95x70-80cm, nonwoven, 4 layers, minimum 42g/m² , waterproof - 2 pcs
 - a. Absorbent zone 85-95x35-45cm
 - b. Adhesive area applied directly to the edge of the field, impregnated in the fabric, without wrinkles, 100% adhesive capacity.
 - c. Latex free, colophony free
3. Celulose absorbent towel 20x25cm - 4pcs
4. Operation bands 5-15x45-55cm, waterproof, absorbent, sticker on all surface, colophony free, for fixation of tubes and instruments - 1 piece
5. Absorbent adhesive drape 165-185x165-185cm, nonwoven, 4 layers, minimum 42g/m² , waterproof -1 piece
 - a. Absorbent zone 60x40cm
 - b. Adhesive area applied directly to the edge of the field, impregnated in the fabric, without wrinkles, 100% adhesive capacity.
 - c. Latex free, colophony free, guiding pictograms
6. Absorbent adhesive drape 230-250x140-160cm, nonwoven, 4 layers, minimum 42g/m², with tube holders, waterproof -1 piece
 - a. Absorbent zone 60x40cm
 - b. Adhesive area applied directly to the edge of the field, impregnated in the fabric, without wrinkles, 100% adhesive capacity.
 - c. Latex free, colophony free, guiding pictograms
7. Absorbent drape for patient 140-160x180-200cm, absorption zone 70-80x185-195cm, 2 layers, waterproof - 1 piece
8. EO sterile