


INJECTORS

OJ 001

IOL INJECTOR FOR B, C CARTRIDGES, TITANIUM




Overall length 160 mm, 6.3 inches

- Screw mechanism
- The lens is advanced through the Cartridge by rotating the plunger
- This results in a controlled defolding of the lens into the capsular bag

OJ 006

IOL INJECTOR FOR D CARTRIDGES, TITANIUM





Overall length 160 mm, 6.3 inches

- Screw mechanism
- The lens is advanced through the Cartridge by rotating the plunger
- This results in a controlled defolding of the lens into the capsular bag

OJ 004

IOL INJECTOR FOR B, C CARTRIDGES, TITANIUM




Overall length 200 mm, 7.8 inches

- Plunger Mechanism with reverse inner spring for easy and efficient one-handed technique implantation
- With ring
- Can be used with IOL loading forceps OF 196

OJ 005

IOL INJECTOR FOR D CARTRIDGE, TITANIUM




Overall length 203 mm, 7.9 inches

- Plunger Mechanism with reverse inner spring for easy and efficient one-handed technique implantation
- With ring
- Can be used with IOL loading forceps OF 196

OJ 007

INJECTOR FOR CAPSULAR TENSION RING



Overall length 161 mm, 6.3 inches

- Stainless hook 0.3 mm dia, 0.85 mm height
- Titanium trunk