





Arterial



Venous

Arterial Tubing(red)

NO.	Description	Qty	Mark
1	Cap for Dialyzer Connector	1	
2	Dialyzer Connector(red)	1	
3	Cap for Female Luer Lock	2	
	Female Luer Lock $\phi 4.1$		
	Clamp(red small)		A
	Branch Tubing: 2.2*4.1*100mm		
4	Cover for Drip Chamber	1	
	Dripper		
5	Pump Segment Connector $\phi 4.1$	1	
	Pump Segment Connector $\phi 5.5$		
	Pump Segment: 8*12*350mm		D
6	Clamp(red small)	1	
	Female Luer Lock $\phi 2.5$		
	Cap for Female Luer Lock		B
	Heparin Line: 1.0*2.5*600mm		C
	Branch Tubing: 2.2*4.1*50mm		
7	Clamp(red small)	1	
	Cap for Female Luer Lock		
	Transducer Protector		E
	Branch Tubing: 3.5*5.5*500mm		
8	Access Port(red)	1	
9	Clamp(red big)	1	
10	Patient Connector	1	
11	Suspended Spike	1	
12	Main Tubing: 4.5*6.8*750mm	1	
13	Main Tubing: 4.5*6.8*600mm	1	
14	Main Tubing: 4.5*6.8*200mm	1	
15	Main Tubing: 4.5*6.8*1600mm	1	

Venous Tubing(blue)

NO.	Description	Qty	Mark
1	Cap for Dialyzer Connector	1	
2	Dialyzer Connector(blue)	1	
3	Cap for Female Luer Lock	1	
	Transducer Protector		
	Clamp(blue small)(blue)		B
	Branch Tubing: 3.5*5.5*500mm		
4	Cover for Drip Chamber	1	
	Drip Chamber (OD:21.5-30mm,LD:130mm)		
	Filter		
5	Access Port(blue)	1	
6	Cap for Female Luer Lock	2	
	Female Luer Lock $\phi 4.1$		
	Clamp(blue small)(blue)		A
	Branch Tubing: 2.2*4.1*100mm		
7	Clamp(blue big)	1	
8	Patient Connector	1	
9	Drainage bag	1	
10	Main Tubing: 4.5*6.8*400mm	1	
11	Main Tubing: 4.5*6.8*500mm	1	
12	Main Tubing: 4.5*6.8*1800mm	1	
	Recirculating Connector	1	

					Prepared by :	Reviewed by :	Description:
					Name		
					Date		
					Approved by QA:		Tubing Sets for Hemodialysis
					Name		
					Date		
3	Rev.	Version:A.2Q	2022/9/5	lai	Fresenius 2008-4008S		Internal ID:
2	Rev.	Version:A.1Q	2022/7/27	lai			BLU006G
1	New Release	Version:AQ	2020/10/12	liu			Blood pathway min. vol. 148.50mL
Rev	Remarks	Version No.	Date	Name			
CARTON SIZE		530*381*249mm		QTY/CTN	24sets		Sterilization R DEHP