



ARES

HIGH PRESSURE SYRINGE

Our quality
is our pride




Interventional Radiology

ARES

Distinguished Service

REMARKABLE CUSTOMER SERVICE

Everything starts with the customer. Treating our customers as our precious business partners, we are dedicated to offer "the best". It is our company principle to deliver the highest level of customer service and assure absolute customer satisfaction. Happy customers are our reward for excellent service.




ARES IS WHERE YOUR QUESTIONS ARE ANSWERED

With our dedicated team well-trained in medical business, we are ready and eager to answer all your questions and provide the information you need. Offering an accessible customer care, Ares is where you can find answers to your questions.




ALWAYS ON TIME

Timeliness is the best in all matters. Ares Medikal embraces the significance of timeliness in delivering customer orders and accepts it as an integral part of customer satisfaction. During two decades in the market, we have developed remarkable skills to meet high demands while assuring the high quality.



OUR COMMITMENT IS HIGH QUALITY

Focusing on research and development and being committed to continuous improvement, we deliver the highest quality products to meet all procedural needs and withstand all procedural challenges.



KNOWLEDGE IS POWER

Our customer-oriented sales team is specialized to provide the clinical knowledge that our customers need. Being in close relation with specialists and interventional radiologists, we will support you to increase your efficiency in the field.

Perfect Syringes for Perfect Contrast Media Injection

ARES

for Tomography



CT-200-A2/DUAL

Injector Model

Medrad Stellant

Volume

200 mL / 200 mL

Content

2 Syringes, 2-Way 150 cm Extension Line W/Check Valves, 2 J tubes, 2 Spikes

Pressure Limit

1200 PSI /84 bar

Material

Polypropylene, Polycarbonate, Latex Free Synthetic Rubber

REF Code: AS8001

- 2 pcs. syringe of 200 ml with transparent material, resistant to pressure of at least 400 psi and FluiDot type markings
- 1 pc Spiral patient connector tube in T with the following characteristics: length 150cm (+-10cm)
- 1 pc Quick fill tube, allows connection to the injector syringes via the T-connector.
- 2 pcs Spikes for