

CAISLEY INTERNATIONAL GMBH – HARDERHOOK 31 – D-46395 BOCHOLT
I.P. Centrul Republican de Diagnostică
Veterinară
str. Murelor, 3,
MD-2051, mun.Chişinău, mun.Chişinău,
Republica Moldova

Ihr Zeichen:
Your Ref:

Unser Zeichen:
Our Ref:
CN/IW

Ihre Nachricht vom
Your message of

Datum
Date
25.11.2023

Compatibility of Universal Applicator

We, CAISLEY International GmbH, hereby declare that the offered eartag models FlexoPlus G/D, MultiFlex L/L and MultiFlex P/P are compatible to the applicator, specified in chapter III of annex 1 and chapter V of annex 2. CAISLEY International GmbH is the manufacturer of this “universal applicator for multiflex ear tags” and we guarantee full compatibility.

The special design of the locking system of the MultiFlex and FlexoPlus causes a clean cut of the ear tissue. The pin of the father part has a circumferential, very hard cutting edge, which together with the mother part cuts the tissue cleanly and does not crush it by pressure. Together with the parallel movement of the CAISLEY Universal Applicator during closure this ensures that no skin is pressed into the mother part, promotes good wound healing and prevents necrosis etc..

Cornelia Nehls
Managing Director


Universal-Ohrmarkenzange

1 Zange für alle Caisley Ohrmarken

Universal Applicator

One Applicator for all Caisley Tags

- FlexoPlus Geno, FlexoPlus
- MultiFlex, FlexoTronic
- PrimaFlex, Duo



Patentierte Parallelschließung

Patented Parallel Closure

Universell einsetzbar
Extrem robust

*Universal Use
Extremely Sturdy*

Produkt Spezifikation

Product specification

Material
Zangenkörper
Zangenköpfe

*Material
Applicator Body
Applicator Heads*

Rostfrei / Rustproof
verzinkter Stahl / Galvanized Steel
GFK / fibre reinforced plastic

Abmessungen
Gewicht

*Dimensions
Weight*

9 x 24 x 3 cm
570 g.

Wechseleinsätze

Changing Inserts

nach Wahl fertig montiert
inserts ready installed as chosen

Grau / Grey




für FlexoPlus-, FlexoPlus Geno-, FlexoTronic-
und Multiflex-Ohrmarken
*for FlexoPlus-, FlexoPlus Geno-, FlexoTronic-
and Multiflex Ear Tags*

Stützringe / Thrust Ring

Grau passt zu grauer und grüner Ohrmarkenspitze
Rot passt zu roter Ohrmarkenspitze
Grey fits to Grey and Green male tip
Red fits to red male tip

Schwarz / Black



Schwarz für PrimaFlex-Ohrmarken
Black for PrimaFlex Ear Tags
Grün für Duo-Ohrmarken
Green for Duo Ear Tags

Grün / Green

VETKIT



Probennahme ohne Tierkennzeichnung
Tissue Sampling without animal identification

Blau / Blue



für Ohrmarken von anderen Herstellern
for Ear Tags from different suppliers