

Shanghai HBM Healthcares, Inc.

#1508-1509,No.211, Golden Port Road, Pudong, Shanghai, China. Tel: 0086-21-58445785 www.hbmchina.com

Medispo[®] Medical gowns

1. Description

- --- designed with classic neck and waist ties, for easily adjustable.
- --- Elastic cuffs.
- ---Color: Blue ,Green and other colors also available
- ---Gowns up to size 3XL, regular and extra long lengths.

2. Composition

100% Spunbounded

3. Specification

Size	REF	Fabric Weight	Length	Packing
S	9630S	35g	110cm	
М	9631M	35g	110cm	
L	9632L	35g	115cm	10Pcs /polybag
XL	9633XL	35g	120cm	100pcs /Carton
XXL	9634XXL	35g	130cm	
XXXL	9635XXXL	35g	140cm	



*Others fabric weight also available.

4. Usage instruction

- ----It is used for medical staff to avoid being contaminated by blood, body fluids, and other infectious substances during contact.
- ---Isolation gown is a two-way isolation to prevent medical staff from being infected or contaminated as well as patients from being infected.
- ---When entering a strictly isolated ward, you need to wear isolation gowns.
- ----Examination and nursing require special isolation of patients. When the work clothes may be contaminated by secretions, excrement, blood or body fluids, isolation clothes shall be worn Medical staff entering the patient wards of infectious diseases that are likely to cause dissemination in the hospital and patients who require special isolation (such as extensive burns, organ transplants, and premature babies) are required to wear isolation gowns.
- ---Single use only.

5. Sterilization Non-Sterile