

Wash Out Cannulas / Viscoelastic Cannulas

Wash Out Cannulas

Code	Gauge and size	Description	Tip
M9280	 19G x 1" 1.1 x 25mm	Washout Cannula. Angled 32°, 8mm from tip (Bishop Harmon). Tip flattened in the horizontal plane	
M9281	 16G x 1" 1.6 x 25mm	Washout Cannula. Angled 30°, 6.5mm from the tip. Tip flattened in the horizontal plane	

Viscoelastic Cannulas

For injection of viscoelastic products into the anterior or posterior chamber.

Code	Gauge and size	Description	Tip
M3926A	 27G x 7/8" 0.4 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M3927	 27G x 7/8" 0.4 x 22mm	Viscoelastic Cannula. Angled at 25°, 8mm from tip	
M2272	 26G x 7/8" 0.45 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M2639A	 25G x 7/8" 0.5 x 22mm	Viscoelastic Cannula. Angled at 35°, 9mm from tip	
M3909	 23G x 7/8" 0.64 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip	
M3910A	 21G x 7/8" 0.81 x 22mm	Viscoelastic Cannula. Angled at 45°, 9mm from tip	
 M3910	 21G x 7/8" 0.81 x 22mm	Viscoelastic Cannula. Straight	
 M3917	 27G x 7/8" 0.4 x 22mm	Viscoelastic Cannula. Angled and formed 11mm from tip	