

Combi stopper



Product Information:

- precisely and safely closes standard ports
- Luer-Lock stopper has only one male tip
- combi stopper with male and female tip especially dedicated
- for securing Luer Lock fitting in infusion and transfusion sets
- single-use
- non-pyrogenic
- latex-free, phthalate-free
- EO sterilized
- packaging: 1 pc./paper-foil

indication

- as a stopper of intravenous cannulas, ports and catheters
- intended to be used in combination with other devices during infusion therapy, protecting against the outflow of blood through the cannula and protecting against contamination and infection



Suzhou Health Plastic Products Co., Ltd
No.699 Ping'an Road Linhu Town Wuzhong
District 215105 Suzhou jiangsu China



MedPath GmbH
Mies-van-der-Rohe-Strasse 8, 80807
Munich, Germany