

**Medispo**<sup>®</sup>  
Premium

We do nothing or, Do Better



## WHAT IS **HBM**

**H HONEST**  
**B BETTER**  
**M MUTUAL**

In HBM-We respect the HONESTY,  
do BETTER, and share MUTUAL benefits.

**Medispo**<sup>®</sup>

BEST IN MEDICAL DISPOSABLE



WE MAKE THE BEST WE CARE THE MOST

 [www.hbmchina.com](http://www.hbmchina.com)

 [gaochen@hbmchina.com](mailto:gaochen@hbmchina.com)

 1508-09, No.211, GoldenPort Road, Shanghai, China




**HBM GROUP**

# Medispo®

## Disposable Face Mask

with ear-loop



Three Layers Design

- Disposable
- Odorless
- 3 layers design
- Comfort and fits perfectly
- Free of fiber glass
- Hypoallergic
- With nose wire support
- Low resistance for better breathing

### PACKING INFORMATION

Packing:	50pc/box, 40box/carton
Net Weight:	7.3kgs/ctn
Gross Weight:	7.9kgs/ctn
Carton Size:	52x40x30cm

### INSTRUCTIONS

1. Put on the mask, hang on the loops to right and left ear.
2. Using hand to mould the nose piece to the shape of your nose, ensure a close fit and good seal.
3. Fit check before using, mask sure no air leaks at the mask edge.



## DISPOSABLE FACE MASK

### 50PCS/BOX

Size:  
 Length 175mm  
 Width folded 95mm  
 Width unfolded 170mm  
 Length of Earloop 170mm

Product Code:  
 TYPE I: 966-EL  
 TYPE II: 966-EL II  
 TYPE IIR: 966-EL IIR

Type I: BFE ≥ 95%  
 Type II: BFE ≥ 98%  
 Type IIR: BFE ≥ 98%

<Bacterial Filtration Efficiency>



# Medispo®

Anti-Fog

## Disposable Medical Face Mask Type II

with Ties



### Type II : BFE ≥ 98%

<Bacterial Filtration Efficiency>




### Three Layers Design

**Layer 1**  
the colored layer is your first line of defense

**Layer 2**  
protects against fluid while maintaining breathability

**Layer 3**  
soft, non-abrasive inner layer makes the mask comfortable on your skin



#### Fluid Resistance

Tests penetration resistance to synthetic blood at three protection levels—120 mm Hg .



#### Bacterial Filtration Efficiency (BFE)

Tests filtration ability against an aerosol containing bacteria 3.0 microns in size.



#### Differential Pressure (ΔP)

Gauges air flow from one side of a mask to the other, assessing breathability.

- Made of non-paper, non-woven material
- Odorless
- Hypoallergenic
- Head straps, of two layers with the same size for perfect fitting on the face
- Ultrasonic Bonding of the mask with the head straps
- Tear resistant.
- Plastic/Wire Nose clip 12cm, for perfect and stable application of the mask on the nose
- Respiratory Resistance - Differential Pressure (Delta P ΔP): <60 (Pa / cm<sup>2</sup>).

#### Dimensions:

Length 175mm  
Width folded 95mm  
Width unfolded 180mm  
controlled-partially lint free



Anti-Fog

### DISPOSABLE MEDICAL FACE MASK TYPE II WITH TIES

50CS/BOX Product Code: 966-TO-AF

#### PACKING INFORMATION

Packing:	50pc/box, 40box/carton
Net Weight:	8.0kgs/ctn
Gross Weight:	9.0kgs/ctn
Carton Size:	52x39x52cm