# **Oral Care Range**

Intersurgical offer a large range of medical toothbrushes with soft moulded bristles and cleaning implements for both adult and paediatric patients. These have been designed with small, patient friendly brush heads which allow the user to gently and effectively reach all areas of the oral cavities.

The OroCare suction products are ready for attachment to suction tubing straight from the packaging without the need for separate handles allowing for quick and efficient usage.











# Foam swabs should only be used to moisten lips



## Foam swabs



78 MHRA reports concerning risk of detachment <sup>1</sup>



Do NOT remove plaque from teeth as effectively as brushing <sup>2</sup>



Should only be used to moisten lips.

## **OroCare™ Aspire suction toothbrush**

With built in direct attachment for the suction line this toothbrush eliminates the need for an additional yankauer. The long slim neck of this brush helps to reach all areas of the oral cavity without interfering with ET tubes and other equipment.



#### OroCare™ Sensitive oral suction wand

This unique suction wand is used for cleaning delicate gums and tissues. The integral moulded soft head provides a safer solution to the commonly used foam swab. The wand allows application of antibacterial solutions whilst gently cleaning the teeth and soft pallet.



#### OroCare™ Mini

Ideal for patients with a challenging oropharyngeal situation. The long slim neck of this brush assists the user in reaching all areas of the oral cavities.





 Otukoya, R. and Doshi, M., 2018. Selecting the right tools for mouth care delivery in hospitals. Nurs Times, 11, pp.18-21.

 Marino PJ et al (2016) Comparison of foam swabs and toothbrushes as oral hygiene interventions in mechanically ventilated patients: a randomised split mouth study.
BMJ Open Respiratory Research; 3: 1, e000150.

