CNWTC

Chongqing New World Trading Co.,Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China Tel: 0086-023-67728212 Fax:0086-023-67728109

<u>Autoclavable Single Channel High Temperature-Resistant Accurate</u> Fixed Adjustable Pipette

Single channel Pipette, also called pipette gun, is a measuring tool that can transfer liquid from the original container to another container within a certain range. It is widely used in biology, chemistry, Medical diagnosis and Laboratory testing etc.

Product Description:

Single channel Pipette, also called pipette gun, is a measuring tool that can transfer liquid from the original container to another container within a certain range. It is widely used in biology, chemistry, Medical diagnosis and Laboratory testing etc.

Product Features:

[Volume]: This Pipette has two types, Fixed or Adjustable volume, a variety of ranges can be selected, you can choose according to your application.

Fixed volume: 1ul,2ul,5ul,10ul,20ul,30ul,50ul,100ul,200ul,500ul,1000ul,1000ul,etc

Adjustable volume: 0.1-2.5ul,0.5-10ul,5-50ul,50-200ul,200-1000ul, 1000-5000ul, 2-10ml etc

[Tolerable temperature]: Resistant to high temperature, the whole product can be autoclaved and high temperature disinfection.

[Lightweight and Ergonomic]: Ergonomic design provides excellent operating experience

Lightweight handle design, Large display window allows for easy volume identification [Highly accuracy]: Every pipette is tested to ensure precision and accuracy. If the volume during use is not accurate, you can re-calibrate.

[Simple maintenance]: Easy calibration and maintenance

Product Specification:

r routet opecification.								
Item	Volume	Type	Packing unit					
6015001	1 ul							
6015002	2 ul							
6015003	2.5 ul							
6015004	5 ul							
6015005	7 ul							
6015006	10 ul							
6015007	15 ul							
6015008	20 ul							
6015009	25 ul							
6015010	30ul							
6015011	40ul							
6015012	50ul							
6015013	75 ul	Fixed	Independent packing					
6015014	80 ul	Tixou	independent packing					
6015015	100 ul							
6015016	200 ul							
6015017	250 ul							
6015018	500 ul							
6015019	1000 ul							



Chongqing New World Trading Co., Ltd

Address: Room 1001, Unit 2, No. 88, Jianxin East Road, Jiangbei District, Chongqing, China Tel: 0086-023-67728212 Fax: 0086-023-67728109

6015020	2000 ul
6015021	2500 ul
6015022	5000 ul
6015023	10000 ul

Item	Volume	Increment	Туре	Packing Unit
	range			
6015024	0.1-2.5 ul	0.05µl		
6015025	0.5-10 ul	0.1µl		
6015026	2-10 ul	0.5µl		
6015027	5-50 ul	0.5µl		
6015028	10-100 ul	1µI		
6015029	20-200ul	1µI	Adjustable	Independent
6015030	50-200 ul	1µl		packing
6015031	100-1000 ul	5µI		
6015032	200-1000 ul	5µI		
6015033	1000-5000 ul	50µl		
6015034	2-10ml	0.1ml		

Product Details:



Packaging Details

One piece in independent packing, then suitable quantity in cartons.

Single package size: 33X13X5 cm Single gross weight: 0.300 kg

Also can customize package as per your request.



Chongqing New World Trading Co., Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China
Tel: 0086-023-67728212 Fax: 0086-023-67728109

Standing pipette rack

Standing pipette rack, can be applied to any specifications of pipette, pipette storage, convenient use of laboratory personnel, can be placed in different environments pipette shelf, single hole, porous, adjustable, quantitative pipette can be very convenient on the vertical pipette shelf.

Introduction:

Standing pipette rack, can be applied to any specifications of pipette, pipette storage, convenient use of laboratory personnel, can be placed in different environments pipette shelf, single hole, porous, adjustable, quantitative pipette can be very convenient on the vertical pipette shelf.

Feature:

Color : Withe
 Material: Plastic

3. Packing: 48*42*45cm, 20pcs/ctn

4. Application: It is used for storing pipettes in laboratory and convenient for placing

pipettes.

Name	Material	Color	Packing	Carton Size (cm)	Weight (kg)
Pipette rack	Plastic	White,blue	20pcs/ctn	48*42*45cm	14.7KG

Detail photos:



CNWTC

Chongqing New World Trading Co., Ltd

Address: Room 1001,Unit 2,No.88,Jianxin East Road, Jiangbei District,Chongqing, China Tel: 0086-023-67728212 Fax:0086-023-67728109

Advantage:

- 1. Pipette bracket adopts strict technological process, suitable for a variety of work occasions, its use can effectively stable, convenient use of pipette, can be suspended at the same time and multi-channel pipette.;
- 2. Pipette bracket has an anti-skid stable base for easy cleaning;
- 3. Can be mounted almost all brands of pipettes;
- 4. Upper center with a handle, easy to move or turn the frame;
- 5. The pipette is permanently attached to the surface of the support and is made of polycarbonate.
- 6. The pipette bracket can be placed with multiple single-channel pipettes (below 1000ul specifications). There are two discharge brackets, the front row bracket is lower than the back row bracket, and the front and back two discharge brackets are placed in each other's position, so that the user can clearly see all the pipettes placed on the bracket.