Medical electric drill RJ-PD (ND-2011)

IMAGE	Technical data
O RUIJIN	Medical electric canulated drill RJ-PD ND-2011
	Working voltage: 14.4V
	Adjustable speed: 0-1000rpm
	Noise: ≤45dB; Drill chuck detachable.
	Output power: 4N.m Output Voltage:100W
	Sterilization method: autoclavable 135° C
	Handpiece weight: 1420g, with battery 1840g
	Battery charger: PJCG001
	Input voltage: AC110V-60Hz and 220V,50HZ
	Output voltage: 14.4V
	Battery charger status: light indicator.
	Charging time: 2-3 hours.
	It can charge one battery each time.
	Re-chargeable battery PJDC1440-2psc
	Voltage: 14.4V 1800mAh
	Fully charged working time: 20 min non-stop
	Aseptic battery transfer kit PJXD001

Address/Tel/Fax: No.33, EastWanchun Road, Wuhu Economic & Technology Development

Zone, Wuhu Anhui China;

Tel: 86-553-2672510 Fax: 86-553-2672511 Website: **www.ruijinmedical.com**



	Key PJYS001-1pcs
O MAN	Sterilization case PJDB006-1 psc
	Drill bits: 5 pcs RJ-RP-TDB
	K wire adapter RJ-RP-KWA1- 1 psc

Date: June 14, 2024

Signature:

Address/Tel/Fax: No.33, EastWanchun Road, Wuhu Economic & Technology Development

Zone, Wuhu Anhui China;

Tel: 86-553-2672510 Fax: 86-553-2672511 Website: www.ruijinmedical.com