

## **CHIRASORB** braided

CHIRASORB braided is a mid-term absorbable, coated, synthetic, multifilament suture made of 90% glycolide and 10% L-lactide. The substances contained in the coating and suture are nonantigenic and noncollagenous.

Sterile Chirasorb braided is available in violet or natural. The product elicits a minimal acute inflammatory reaction in tissues.

Progressive loss of tensile strength and absorption of Chirasorb braided occurs by hydrolysis, when the polymer is degraded to glycolic and lactic acids, which are subsequently absorbed and metabolised by the body.



Indication: General soft tissue approximation and/or ligation in all specialties, including use in

Ophthalmic surgery.

Contraindication: Chirasorb braided should not be used in tissues subject to expansion, stretching or

distension, or which might require long-term mechanical support.

**Sizes:** 10/0 to 6 USP

Structure: Multifilament

Absorption method: Hydrolysis

**Strength Retention** 

Profile: MID-TERM HEALING (35 days)

1 Day 100% 14 Days 75% 21 Days 50%

**Mass Absorption** 

profile: 56 - 70 days

UK07

65



## CHIRASORB braided

U	n	U	1
	7	2	1

_		

Needle	EP	USP	Length	Colour	Cat. No	Other information				
HRT: ½ circle, round boo	dy needle - ta	percut								
HRT 17	1.5	4/0	1 x 45 cm	Violet	LV 0566					
	1.5 2	4/0 3/0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0561 LV 0562					
HRT 20	1	5/0	1 x 45 cm	Violet	LV 0568					
HRT 22	2	3/0	1 x 75 cm	Violet	LV 0572					
HRT 25	2	3/0	1 x 75 cm	Violet	LV 0581					
HRT 31	3	2/0	1 x 75 cm	Violet	LV 0606					
(HRT 37) \	<mark>3</mark> 3.5	<mark>2/0</mark> 0	1 x 75 cm 1 x 75 cm	Violet Violet	LV 0592 LV 0593					
	3 3.5 4	2/0 0 1	1 x 90 cm 1 x 90 cm 1 x 90 cm	Violet Violet Violet	LV 0597 LV 0598 LV 0599					
HRT 40	3.5	0	1 x 75 cm	Violet	LV 0653					
HRT 45	5	2	1 x 75 cm	Violet	LV 0665					
HRN: ½ circle, round body needle - blunt point										
HRN 50	4 5	1 2	1 x 90 cm 1 x 90 cm	Violet Violet	LV 0686 LV 0675					