



INJECTION NEEDLE


WHERE ERGONOMICS MEET RELIABILITY


MICRO-TECH's disposable injection needle enhances clinical practice due to its ergonomic handle and improved needle sharpness. The automatic "pen-style" handle allows single-handed use. For even more precise punctures, the needle has varying stiffness levels that can be gradually adjusted. When the needle is exposed, the quick-release

button snaps the needle back into the catheter to avoid damage to the working channel. The patented ratchet mechanism provides audible confirmation of a drawn needle, providing an effective exchange between the physician and nurse.

SPECIFIC CHARACTERISTICS

- Ergonomic handle with one-handed operation
- Patented ratchet mechanism with audible confirmation
- Adjustable stiffness levels
- Quick-release button for a safe needle retraction
- Enhanced needle sharpness for precise penetration
- Different needle diameters and lengths
- Easy attachment via a secure Luer-Lock connection





Needle retracted

Needle extended



QUICK-RELEASE
BUTTON

SPECIFICATIONS

REF	Needle length mm	Needle Ø	Working channel Ø mm	Length mm
INJECTION NEEDLE - POLYPECTOMY/EMR				
IN12-195232302	5	19 G = 1.1 mm	2.8	2300
IN12-225232302	5	22 G = 0.7 mm	2.8	2300
IN12-255232302	5	25 G = 0.5 mm	2.8	2300
INJECTION NEEDLE - ESD				
IN12-224232302	4	22 G = 0.7 mm	2.8	2300
IN12-254232302	4	25 G = 0.5 mm	2.8	2300

Packaging unit: 10 pieces