



GIMA

LARYNGOSCOPE SET 4 MC INT. BL. 1-2-3-4 - adult

Code:	34303
Category:	Standard laryngoscope sets
Unit of sale:	1 pc.
Minimum order:	1
Type:	Medical device
Class:	I
NSIS:	2535677
CND:	Z12021003
EAN13:	8023279343038



Description: STANDARD RECHARGEABLE SETS (HANDLE + BLADES)
• Adult set (4 blades Mc INTOSH 1-2-3-4)

Great features of the unique GIMA blades and handles

- standard 2.5 V vacuum bulb
- light ergonomic, satin finish handle to reduce glare
- Stainless Steel contacts guarantee long working life
- easy sterilization ETO or steam
- full range of Mc Intosh & Miller blades, latex free
- ISO 7376 fittings allow use of blades with existing handles
- set is packaged in an anti-shock case and in cardboard box
- 2.5 V handles work with both alkaline batteries or rechargeable batteries
- good light transmission of 3,500 LUX with 2.5 V handles
- autoclavable up to 134° for 5 minutes (approx. 2,000 times)
- non magnetic blades

DECLARAȚIE DE CONFORMITATE

Societatea GIMA S.P.A., cu sediul administrativ în Gessate (MI), Via Marconi 1 și sediul social în Milano, Via Tommaso Grossi 2, în calitate de producător al dispozitivului medical:

Număr unic de înregistrare GIMA (SRN):

Dispozitiv medical (nume și denumire)	Cod	Cod UDI-DI de bază
SET LARINGOSCOPIE 4 LAME - Lame Mc Intosh 1-2-3-4	34303	8023279Z120210030000000 MQ

Clasă de risc I (nesteril), conform reglementării 13 Anexa VIII la Regulamentul (UE) 745/2017 (MDR), declară, pe propria răspundere exclusivă, că acest dispozitiv:

- a fost realizat în conformitate cu cerințele esențiale și cu prevederile Regulamentului (UE) 745/2017 (MDR), așa cum reiese din dosarul tehnic arhivat la sediul firmei;
- nu s-au utilizat specificații comune pentru conformitatea dispozitivului medical sus-menționat;
- a fost realizat în conformitate cu directiva 2011/65/UE (cu modif. și completările ulterioare) privind interzicerea utilizării anumitor substanțe periculoase în echipamentele electrice și electronice.

Gessate, 5/28/2021

GIMA S.p.A.
Reprezentant
legal(Nicola



Manzoni)