



Bone Cement

With Antibiotic

Radiopaque Sterile
Surgical Acrylic

41G



A pouch containing an envelope with 41 g of radiopaque

Gentamicin-loaded powder sterilized by ethylene oxide

Polymethyl Methacrylate 35,58 g (86,78%)

Barium Sulphate Ph. Eur. 4 g (9,76%)

Benzoyl Peroxide 0,42 g (1,02%)

Gentamicin sulfate 1 g (2,44%)

An ampoule with 20 ml of liquid sterilized by ultrafiltration


Methyl Methacrylate 19,6 ml (98%)

N,N dimethyl p-toluidine 0,4 ml(2%)

Hydroquinone 18-20 ppm



Alcasis Genta Working Duration



Mpa Flexural Strength

