

Alcohol Swab Dispenser QUICKPAD®

Practical and efficient

The QUICKPAD® alcohol swab dispenser, being sterile and physiologically verified as harmless, is ideal for cleaning and disinfecting the skin. The active agent 2-propanol acts as an effective, "mild on skin" disinfecting agent. The patented "swab tear-off" system allows for the single use of the swab. This makes QUICKPAD® not only economical but also efficient in its use.

Quality and Endurance

The lid of the QUICKPAD® alcohol swab dispenser container seals air-tight, keeping the alcohol swabs moist and sterile. As a result, the swab dispenser has a particularly long shelf life of 24 months.

Please note that for pharmaceutical QUICKPAD®, the corresponding safety instruction is applicable



Characteristics

- Single swab tear-off system allows for an efficient and economical use;
- Fill Level control due to transparent container;
- Retains moisture and sterility;
- Physiologically harmless non-woven swab;
- 150 swabs (44 x 44 mm);
- Suitable for self-application by the patient;
- Ready for use;
- Made in Germany Quality, CE Label;

Range of Application

QUICKPAD® is supplied ready for use and with its simple operation, can be used by specialists as well as patients for disinfecting the skin. An ideal product for diabetics and other users of subcutaneous self-injection syringes.