

# Certificate of Analysis



## Prep-n-Go™ Buffer

Product No. **4471406**  
Lot No. **2407053**  
Date of Manufacture **15AUG2024**  
Expiration Date **02JAN2026**

### KIT CONFIGURATIONS:

PART NO.	COMPONENT DESCRIPTION	LOT NO.
4471434	(1) Bottle Prep-n-Go™ Buffer	2407050

### SPECIFICATIONS:

### RESULT

QC Method for Prep-n-Go™ Buffer: The analytical QC results are obtained from PN 4467081 BULK, Prep-n-Go™ Buffer. The human DNA contamination results are obtained from P/N 4471434 BOTTLE, Prep-n-Go™ Buffer and are analyzed on an Applied Biosystems Genetic Analyzer using the AmpFℓSTR® Identifiler® Direct PCR Amplification Kit.

No human STR profile detected using AmpFℓSTR® Identifiler® Direct PCR Amplification Kit at 28 cycles

Pass

pH (8.15 - 8.30)

Pass

Conductivity (1.45 - 1.88 mS/cm)

Pass

Refractive Index (1.3330 - 1.3340)

Pass

***For Research Use Only. Not for use in diagnostic procedures.***

Manufactured in compliance with our ISO 13485 certified quality management system.  
Site: Warrington, UK

Life Technologies  
5781 Van Allen Way  
Carlsbad, CA, USA 92008  
[www.thermofisher.com](http://www.thermofisher.com)  
For inquiries, contact us at

[cofarequests@thermofisher.com](mailto:cofarequests@thermofisher.com)

**Quality Assurance**  
**Issued 19AUG2024**

Doc: 4471046 Rev.E