

# Manifolds



# **PRODUCT FEATURES**

- Compatible with most infusion pumps, syringes and intravenous infusion tubes to meet complex clinical needs
- Excellent pressure resistance and sealing performance
- Leak free and smooth handling
- On / Off types are available for different user habits

# **ONE PORT**

Catalogue No.	Product description/Pressure
BM-MF1-ON	"ON" 500psi
BM-MF1-OFF	"OFF" 500psi

# **TWO PORTS**

Catalogue No.	Product description/Pressure
BM-MF2-ON	"ON" 500psi
BM-MF2-OFF	"OFF" 500psi

#### **THREE PORTS**

Catalogue No.	Product description/Pressure
BM-MF3-ON	"ON" 500psi
BM-MF3-OFF	"OFF" 500psi



# **Radial Artery Tourniquets**



# **INTRODUCTION**

The device is used to achieve hemostasis after the radial arterial catheterization.

#### **PRODUCT FEATURES**

- Unique design for radial artery hemostasis post percutaneous puncture
- Transparent fixed plates ensures clear visual control
- Soft and comfortable adjustable band

#### **PRODUCT COMPOSITION**

The device consists of a fixed plate, adjustable band and elastomer

Catalogue No.	Product description
BM-RAT1	Cord button
BM-RAT2	Rotary

10 www.brosmed.com 11