

Vacutainer Holder



Vacutainer Holder

Blood collection needle holder is compatible with Multi-Sample Needle to collecting blood.

Lock the multi-sample needle short end which is with latex cover directly into holder. Insertions finish when you push the white part of the snap.

After collection finished, press the green color button on holder, needle automatically discharge. Collector can easily finish collection without touch of needle and avoid mistake puncture.

It is non-dangerous products, non-flammable, non-explosive, and can be stored at room temperature.

Cat. No.	Description	Qty/Case(pcs)
632201	Vacutainer holder	4000
632202	Vacutainer holder (safety type)	4000