Infusion/Transfusion

Pressure Infusion Cuffs



Pressure Infusion Cuffs

The safe press can be used for liquid bags up to 500 resp. 1.000 ml filling volume. It consists of one cuff with integrated bag holder, one connecting tube with manometer and one pumpball. The cuff is used for pressure-supported infusions and transfusions.

Sterile: No

REF	Capacity	PU
43.018.00.012	500 ml	20 pcs
43.019.00.012	1.000 ml	20 pcs



Pressure Infusion Cuffs, reusable

The safe press can be used for liquid bags up to 500 resp. 1.000 ml filling volume. It consists of one cuff with integrated bag holder, one connecting tube with manometer and one pumpball with outlet valve. The cuff is used for pressure-supported infusions and transfusions. The cuff is reusable.

Allergenic potential: Latex-free

REF	Capacity	PU
43.018.00.013	500 ml	6 pcs
43.019.00.013	1.000 ml	6 pcs

02-03/2019

Product Data Sheet



Pressure Infusion Cuff

REF: 43.018.00.013

Unit: 1 pc

Packing Unit: 6 pcs.

Description

The re-usable Pressure infusion cuff "Safe Press" is suitable for up to 500 ml liquid collection bags. It contains a cuff with an integrated hanger, a connection tube with aneroid gauge and a hand bulb with integrated bottom valve.

The cuff is used for pressure-supported infusions and transfusions.

The cuff is re-usable.

The cuff does not contain latex.

Material Composition

Product: bag: Nylon

hanger: ABS tube: PVC aneroid gauge: ABS valve: PVC bottom valve: PVC

Packing: individual packing: PE-Bag

outer case: corrugated cardboard

Sterilisation

Non-sterile

