

## Flow Regulator | Code : 1010

### •Description

• Flow regulator is used to administer Intravenous fluid and medicines into human circulating system.

### •Features

1- Precise flow rate control of IV fluids with range of 5 to 250 mL/hr.



## Three Way Stop Cock | Code : 1006

### •Description

- The Three Way Stop Cock family is designed to provide access into I.V. Cannula system for administration of two fluids or drugs.
- For pressure up to 4.5 bar (65 psi).
- 360-degree rotation.

### •Features

1- Rotating male luer and two female luer connectors.  
2- Arrow indication on the handle to indicate the direction of flow.





### Three Way Stop Cock with Low Pressure Extension Tube

#### Description

Code : 1008-LP

- The Three Way Stop Cock family is designed to provide access into I.V. Cannula system for administration of two fluids or drugs at a time, based on the application and requirements.
- For pressure up to 4.5 bar (65 psi).
- 360-degree rotation.

#### Features

- 1- Rotating male luer and two female luer connectors.
- 2- Arrow indication on the handle to indicate the direction of flow.
- 3- Tube Length:  
(7cm - 300 cm).