## CINCINNATI SUB-ZERO Normothermia System

# Gelli-Roll®

The revolutionary and simple-to-use Gelli-Roll® combines patient temperature management and comfort. The Gelli-Roll® is a reusable water blanket encapsulated in polymer. It may be used before, during, or after surgery to effectively keep your patient normothermic.



- Easy to clean with disinfectants
- Reduces set-up and breakdown time
- Supports most patients, infant to adult
- Allows better access to the patient then blankets





### Gelli-Roll®

#### Reusable Hypo-Hyperthermia Blanket

Ordering Information		
Cat. No.	Description	Quanity
195P	Gelli-Roll® - Adult 73.5" x 21.5" x 0.6" (186.7cm x 54.6cm x 1.6cm)	ea.
194P	Gelli-Roll® - Pediatric 31.3" x 23.7" x 0.6" (79.4cm x 60.3cm x 1.6cm)	ea.
193P	Gelli-Roll® - Infant 20" x 13" x 0.6" (50.8cm x 33cm x 1.6cm)	ea.
Hyper-Hypothermia units		
233	Blanketrol® III	115v, 230v
2225	Blanketrol®II	115v
111W	Norm-O-Temp®	230v, 240v
Reusable Connecting Hose (must purchase 286 hose to use Gelli-Roll with hyper-hypothermia units)		
286	9' Reusable Connecting Hose (2.7m)	ea.
286-18	18' Reusable Connecting Hose (5.5m)	ea.
286-27	27' Reusable Connecting Hose (8.2m)	ea.



#### **Product Specifications:**

**Dimensions:** 

Adult - 195P:

73.5in x 21.5in x 0.6in (186.7cm x 54.61cm x 1.6cm)

Weight: 35 lb. (15.9 kg)

Pediatric - 194P:

31.3in x 23.7in x 0.6in (79.38cm x 60.3cm x 1.6cm)

Weight: 25 lb (11.34 kg)

Infant - 193P:

20in x 13in x 0.6in (50.8cm x 33.02cm x 1.6cm)

Weight: 15 lb (6.8 kg)



Phone: 513-772-8810 Toll Free: 800-989-7373 Fax: 513-772-9119

www.cszmedical.com

